

# **SONOACE X8 Ultrasound System**

## **DICOM Conformance Statement**

**Rev 1.0**

Created: March 19, 2007  
Last Updated: May 3, 2007



## **0 COVER PAGE**

Company Name: MEDISON CO., LTD

Produce Name: SONOACE X8

Version: 1.0

Internal Document Number:

Date: May 3, 2007

# 1 CONFORMANCE STATEMENT OVERVIEW

SONOACE X8 implements the necessary DICOM services to download worklists from information systems, save acquired US images and Structured Reports to a network storage device, CD or DVD, print to a networked hardcopy device and inform the information system about the work actually done.

Table 1-1 provides an overview of the network services supported by SONOACE X8.

**Table 1-1  
NETWORK SERVICES**

| <b>SOP Classes</b>                   | <b>User of Service (SCU)</b> | <b>Provider of Service (SCP)</b> |
|--------------------------------------|------------------------------|----------------------------------|
| <b>Transfer</b>                      |                              |                                  |
| Ultrasound Image Storage             | Yes                          | No                               |
| Ultrasound Multi-frame Image Storage | Yes                          | No                               |
| Comprehensive SR                     | Yes                          | No                               |
| <b>Workflow Management</b>           |                              |                                  |
| Modality Worklist                    | Yes                          | No                               |
| Storage Commitment Push Model        | Yes                          | No                               |
| Modality Performed Procedure Step    | Yes                          | No                               |
| <b>Print Management</b>              |                              |                                  |
| Basic Grayscale Print Management     | Yes                          | No                               |
| Basic Color Print Management         | Yes                          | No                               |

Table 1-2 provides an overview of the Media Storage Application Profiles supported by SONOACE X8.

**Table 1-2  
MEDIA SERVICES**

| <b>Media Storage Application Profile</b> | <b>Write Files (FSC or FSU)</b> | <b>Read Files (FSR)</b> |
|--|---------------------------------|-------------------------|
| <b>DVD – Rewritable</b>                  |                                 |                         |
| STD-US-SC-MF-DVD                         | Yes                             | No                      |

## 2 TABLE OF CONTENTS

|           |  |    |
|-----------|--|----|
| 0         | COVER PAGE.....  | 2  |
| 1         | CONFORMANCE STATEMENT OVERVIEW.....  | 3  |
| 2         | TABLE OF CONTENTS .....  | 4  |
| 3         | INTRODUCTION.....  | 8  |
| 3.1       | REVISION HISTORY.....  | 8  |
| 3.2       | AUDIENCE .....   | 8  |
| 3.3       | REMARKS.....   | 8  |
| 3.4       | DEFINITIONS, TERMS AND ABBREVIATIONS.....                                  | 9  |
| 3.5       | REFERENCES.....  | 9  |
| 4         | NETWORKING .....   | 10 |
| 4.1       | IMPLEMENTATION MODEL .....   | 10 |
| 4.1.1     | Application Data Flow .....  | 10 |
| 4.1.2     | Functional Definition of AE's.....   | 11 |
| 4.1.2.1   | Functional Definition of Storage Application Entity.....                   | 11 |
| 4.1.2.2   | Functional Definition of Workflow Application Entity.....                  | 11 |
| 4.1.2.3   | Functional Definition of Hardcopy Application Entity .....                 | 12 |
| 4.1.3     | Sequencing of Real-World Activities .....                                  | 13 |
| 4.2       | AE SPECIFICATIONS.....   | 15 |
| 4.2.1     | Storage Application Entity Specification.....                              | 15 |
| 4.2.1.1   | SOP Classes .....  | 15 |
| 4.2.1.2   | Association Policies .....   | 15 |
| 4.2.1.2.1 | General.....   | 15 |
| 4.2.1.2.2 | Number of Associations .....   | 15 |
| 4.2.1.2.3 | Asynchronous Nature.....   | 16 |
| 4.2.1.2.4 | Implementation Identifying Information.....                                | 16 |
| 4.2.1.3   | Association Initiation Policy .....  | 16 |
| 4.2.1.3.1 | Activity – Send Images and Structured Reports and Requests Commitment..... | 16 |
| 4.2.1.4   | Association Acceptance Policy .....  | 23 |
| 4.2.1.4.1 | Activity – Receive Storage Commitment Response.....                        | 23 |
| 4.2.2     | Workflow Application Entity Specification.....                             | 25 |
| 4.2.2.1   | SOP Classes .....  | 25 |
| 4.2.2.2   | Association Establishment Policy .....                                     | 26 |

|                              |  |           |
|------------------------------|--|-----------|
| 4.2.2.2.1                    | General.....   | 26        |
| 4.2.2.2.2                    | Number of Associations.....                                    | 26        |
| 4.2.2.2.3                    | Asynchronous Nature.....                                       | 26        |
| 4.2.2.2.4                    | Implementation Identifying Information.....                    | 26        |
| 4.2.2.3                      | Association Initiation Policy.....                             | 27        |
| 4.2.2.3.1                    | Activity – Worklist Update.....                                | 27        |
| 4.2.2.3.2                    | Activity – Acquire Images.....                                 | 31        |
| 4.2.2.4                      | Association Acceptance Policy.....                             | 36        |
| 4.2.3                        | Hardcopy Application Entity Specification.....                 | 37        |
| 4.2.3.1                      | SOP Classes.....   | 37        |
| 4.2.3.2                      | Association Policies.....                                      | 37        |
| 4.2.3.2.1                    | General.....   | 37        |
| 4.2.3.2.2                    | Number of Association.....                                     | 37        |
| 4.2.3.2.3                    | Asynchronous Nature.....                                       | 37        |
| 4.2.3.2.4                    | Implementation Identifying Information.....                    | 38        |
| 4.2.3.3                      | Association Initiation Policy.....                             | 38        |
| 4.2.3.3.1                    | Activity – Film Images.....                                    | 38        |
| 4.2.3.4                      | Association Acceptance Policy.....                             | 46        |
| 4.3                          | NETWORK INTERFACE.....   | 47        |
| 4.3.1                        | Physical Network Interface.....                                | 47        |
| 4.4                          | CONFIGURATION.....   | 47        |
| 4.4.1                        | AE Title/Presentation Address Mapping.....                     | 47        |
| 4.4.1.1                      | Local AE Titles.....   | 47        |
| 4.4.1.2                      | Remote AE Title/Presentation Address Mapping.....              | 47        |
| 4.4.1.2.1                    | Storage.....   | 47        |
| 4.4.1.2.2                    | Workflow.....  | 48        |
| 4.4.1.2.3                    | Hardcopy.....  | 48        |
| 4.4.2                        | Parameters.....  | 48        |
| <b>5</b>                     | <b>MEDIA INTERCHANGE.....</b>                                  | <b>50</b> |
| 5.1                          | IMPLEMENTATION MODEL.....                                      | 50        |
| 5.1.1                        | Application Data Flow.....                                     | 50        |
| 5.1.2                        | Functional Definition of AEs.....                              | 50        |
| 5.1.2.1                      | Functional Definition of Offline-Media Application Entity..... | 50        |
| 5.1.3                        | Sequencing of Real-World Activities.....                       | 50        |
| 5.1.4                        | File Meta Information Options.....                             | 50        |
| 5.2                          | AE SPECIFICATIONS.....   | 51        |
| SONOACE X8 DICOM Conformance |  | 5         |

|           |  |           |
|-----------|--|-----------|
| 5.2.1     | Offline-Media Application Entity Specification.....        | 51        |
| 5.2.1.1   | File Meta Information for the Application Entity.....      | 51        |
| 5.2.1.2   | Real-World Activitiesdynamix0.....                         | 51        |
| 5.2.1.2.1 | Activity – Export to DVD.....                              | 51        |
| <b>6</b>  | <b>SUPPORT OF CHARACTER SETS.....</b>                      | <b>53</b> |
| <b>7</b>  | <b>SECURITY .....</b>                                      | <b>54</b> |
| <b>8</b>  | <b>ANNEXES.....</b>  | <b>55</b> |
| 8.1       | IOD CONTENTS .....   | 55        |
| 8.1.1     | Created SOP Instances.....                                 | 55        |
| 8.1.1.1   | US or US Multiframe Image IOD .....                        | 55        |
| 8.1.1.2   | Comprehensive Structured Report IOD .....                  | 56        |
| 8.1.1.3   | Common Modules.....  | 56        |
| 8.1.1.4   | US or US Multiframe Image Module.....                      | 60        |
| 8.1.1.5   | Comprehensive Structured Report Modules .....              | 64        |
| 8.1.2     | Used Fields in received IOD by application.....            | 67        |
| 8.1.3     | Attribute mapping.....                                     | 67        |
| 8.1.4     | Coerced/Modified Fields .....                              | 68        |
| 8.2       | DATA DICTIONARY OF PRIVATE ATTRIBUTES.....                 | 69        |
| 8.3       | CODED TERMINOLOGY AND TEMPLATES .....                      | 69        |
| 8.4       | GRAYSCALE IMAGE CONSISTENCY .....                          | 69        |
| 8.5       | STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES..... | 69        |
| 8.5.1     | US OR US MULTIFRAME IMAGE STORAGE SOP CLASS.....           | 69        |
| 8.6       | PRIVATE TRANSFER SYNTAXES.....                             | 69        |
| <b>9</b>  | <b>STRUCTURED REPORT TEMPLATES.....</b>                    | <b>70</b> |
| 9.1       | TEMPLATES USED IN SONOACE X8.....                          | 70        |
| 9.1.1     | OB-GYN STRUCTURED REPORT TEMPLATE .....                    | 71        |
| 9.1.1.1   | OB-GYN Ultrasound Procedure Report (TID 5000) .....        | 71        |
| 9.1.1.2   | Observation Context (TID 1001).....                        | 72        |
| 9.1.1.3   | Patient Characteristics (TID 5001).....                    | 73        |
| 9.1.1.4   | OB-GYN Procedure Summary Section (TID 5002).....           | 75        |
| 9.1.1.5   | Fetal Biometry Ratio Section (TID 5004).....               | 78        |
| 9.1.1.6   | Fetal Biometry Section (TID 5005) .....                    | 79        |
| 9.1.1.7   | Long Bones Section (TID 5006).....                         | 83        |
| 9.1.1.8   | Fetal Cranium Section (TID 5007).....                      | 85        |

|          |  |     |
|----------|--|-----|
| 9.1.1.9  | Fetal Biophysical Profile Section (TID 5009).....                              | 87  |
| 9.1.1.10 | Early Gestation Section (TID 5011).....  | 88  |
| 9.1.1.11 | Amniotic Sac Section (TID 5010).....   | 91  |
| 9.1.1.12 | Pelvis and Uterus Section (TID 5015).....                                      | 92  |
| 9.1.1.13 | Ovaries Section (TID 5012).....  | 94  |
| 9.1.1.14 | Follicles Section - Left (TID 5013).....                                       | 96  |
| 9.1.1.15 | Follicles Section – Right (TID 5013).....                                      | 97  |
| 9.2      | DCMR CONTEXT GROUPS USED IN SONOACE X8.....                                    | 99  |
| 9.2.1    | Gestational Age Equations and Tables (Context Group 12013).....                | 99  |
| 9.2.2    | OB Fetal Body Weight Equations and Tables (Context ID 12014).....              | 100 |
| 9.2.3    | Fetal Growth Equations and Tables (Context ID 12015).....                      | 101 |
| 9.2.4    | Estimated Fetal Weight Percentile Equations and Tables (Context ID 12016)..... | 102 |

## 3 INTRODUCTION

### 3.1 REVISION HISTORY

| Document Version    | Date of Issue  | Author  | Description                   |
|---------------------|----------------|---------|-------------------------------|
| 1.0 Draft           | March 20, 2007 | MEDISON | Draft Version                 |
| 1.0 Draft Check     | May 3, 2007    | MEDISON | Checking Version of the draft |
| 1.0 Trial-Use Draft | May 3, 2007    | MEDISON |                               |

### 3.2 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

### 3.3 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with MEDISON and other vendor's Medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. However, it is not guaranteed to ensure by itself the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between MEDISON and non – MEDISON equipment.
- Test procedures should be defined to validate the desired level of connectivity.
- The DICOM Standard will evolve to meet the users' future requirements. MEDISON is actively involved in developing the standard further and therefore reserves the right to make changes to its products or to discontinue their delivery.

### 3.4 DEFINITIONS, TERMS AND ABBREVIATIONS

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM Standard.

Abbreviations and terms are as follows:

|      |   |
|------|---|
| AE   | DICOM Application Entity                    |
| AET  | Application Entity Title                    |
| ASCE | Association Control Service Element         |
| CD-R | Compact Disk Recordable                     |
| FSC  | File-Set Creator                            |
| FSU  | File-Set Updater                            |
| FSR  | File-Set Reader                             |
| IOD  | (DICOM) Information Object Definition       |
| ISO  | International Standard Organization         |
| MPPS | Modality Performed Procedure Step           |
| MSPS | Modality Scheduled Procedure Step           |
| R    | Required Key Attribute                      |
| O    | Optional Key Attribute                      |
| PDU  | DICOM Protocol Data Unit                    |
| SCU  | DICOM Service Class User (DICOM client)     |
| SCP  | DICOM Service Class Provider (DICOM server) |
| SOP  | DICOM Service-Object Pair                   |
| U    | Unique Key Attribute                        |


### 3.5 REFERENCES

[DICOM] Digital Imaging and Communications in Medicine (DICOM), NEMA PS 3.1–3.18, 2007

## 4 NETWORKING

### 4.1 IMPLEMENTATION MODEL

#### 4.1.1 Application Data Flow



**Figure 4.1-1**  
**APPLICATION DATA FLOW DIAGRAM**

- The Storage Application Entity sends images, Structured Reports and requests Storage Commitment to a remote AE. It is associated with the local real-world activities "Send Images & Reports" and "Request Commitments". Methods to send images depend on user configuration, "Batch", "Send As You Go" or

“Manual”. “Manual” mode is performed upon user request for each study or for specific images selected. “Batch” mode starts to send images at End Exam for each study. “Send As You Go” mode starts when the first image is acquired for each study and images are transferred immediately after acquisition.

Structured Reports are only sent at End Exam for each study.

If the remote AE is configured as an archive device, the Storage AE will request Storage Commitment and if a commitment is successfully obtained, it will record this information in the local database and displayed it in the Exam List.

- The Workflow Application Entity receives Worklist information from and sends MPPS information to a remote AE. It is associated with the local real-world activities “Update Worklist” and “Acquire Images”. When the “Update Worklist” local real-world activity is performed, the Workflow Application Entity queries a remote AE for worklist items and provides the set of worklist items matching the query request. “Update Worklist” is performed as a result of an operator request or can be performed automatically at specific time intervals. When the “Acquire Images” local real-world activity is performed, the Workflow Application Entity creates and updates Modality Performed Procedure Step instances managed by a remote AE. Acquisition of images will result in automated creation of an MPPS Instance. Completion of the MPPS is performed at End Exam for each study.
- The Hardcopy Application Entity prints images on a remote AE (Printer). It is associated with the local real-world activity “Film Images”. Methods to film Images depend on user configuration and are equal to the Sending images’ of the Storage Application Entity.

## **4.1.2 Functional Definition of AE’s**

### **4.1.2.1 Functional Definition of Storage Application Entity**

The existence of a send-job with associated network destination will activate the Storage AE. An association request is sent to the destination AEs and upon successful negotiation of a Presentation Context, the image transfer is started. If the association cannot be opened, the related send-job is set to an error state and can be restarted by the user via DICOM manager interface or automatically. An automatic retry (retry interval, retry count) can be configured using the Setup/DICOM Menu.

### **4.1.2.2 Functional Definition of Workflow Application Entity**

Worklist Update attempts to download a Worklist from a remote node. If the Workflow AE establishes an association to a remote AE, it will transfer all matching worklist items via the open Association. By default,


Worklist Update use "US" for Modality, current date for Scheduled Procedure Step Start Date and blank for Scheduled Station AE-Title as query parameters. The results will be displayed in a separate list, which will be cleared with the next Worklist Update.

The Workflow AE performs the creation of an MPPS Instance automatically whenever the first image is acquired for each study. The MPPS "Complete" states can only be set by "End Exam" for each study.


#### **4.1.2.3 Functional Definition of Hardcopy Application Entity**

The existence of a print-job will activate the Hardcopy AE. An association is established with the printers and the printer's status determined. If the printer is operating normally, the film sheets described within the print-job will be printed. If the printer is not operating normally, the print-job will set to an error state and can be restarted by the user via DICOM manager interface or automatically. An automatic retry (retry interval, retry count) can be configured using the Setup/DICOM Menu.

### 4.1.3 Sequencing of Real-World Activities



**Figure 4.1-2**  
**SEQUENCING CONSTRAINTS – SEND AS YOU GO**



**Figure 4.1-3**  
**SEQUENCING CONSTRAINTS – BATCH MODE**

Under normal scheduled workflow conditions, the sequencing constraints are illustrated in Figure 4.1-2 and Figure 4.1-3.

Other workflow situations (e.g. unscheduled procedure steps) will have other sequencing constraints. Printing could equally take place after the images acquired have been stored. Printing could be omitted completely if no printer is connected or hardcopies are not required.

## 4.2 AE SPECIFICATIONS

### 4.2.1 Storage Application Entity Specification

#### 4.2.1.1 SOP Classes

SONOACE X8 provides Standard Conformance to the following SOP Classes:

**Table 4.2-1  
SOP CLASSES FOR AE STORAGE**

| SOP Classes                             | SOP Class UID                 | SCU | SCP |
|---|-------------------------------|-----|-----|
| Ultrasound Image Storage                | 1.2.840.10008.5.1.4.1.1.6.1   | Yes | No  |
| Ultrasound Multi-frame Image Storage    | 1.2.840.10008.5.1.4.1.1.3.1   | Yes | No  |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.33 | Yes | No  |
| Storage Commitment Push Model           | 1.2.840.10008.1.20.1          | Yes | No  |
| Verification                            | 1.2.840.10008.1.1             | Yes | Yes |

#### 4.2.1.2 Association Policies

##### 4.2.1.2.1 General

The DICOM Standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-2  
DICOM APPLICATION CONTEXT FOR AE STORAGE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 4.2.1.2.2 Number of Associations

SONOACE X8 can initiate one or more Associations at a time for each destination to which a transfer request is being processed in the active job queue list.

**Table 4.2-3**

### NUMBER OF ASSOCIATIONS INITIATED FOR AE STORAGE

|   |           |
|---|-----------|
| Maximum number of simultaneous Associations | Unlimited |
|---|-----------|

SONOACE X8 accepts Associations to receive N-EVENT\_REPORT notifications for the Storage Commitment Push Model SOP Class.

**Table 4.2-4**

### NUMBER OF ASSOCIATIONS ACCEPTED FOR AE STORAGE

|   |           |
|---|-----------|
| Maximum number of simultaneous Associations | Unlimited |
|---|-----------|

#### 4.2.1.2.3 Asynchronous Nature

SONOACE X8 does not support asynchronous communications (multiple outstanding transactions over a single Association).

**Table 4.2-5**

### ASYNCHRONOUS NATURE AS A SCU FOR AE STORAGE

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

#### 4.2.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-6**

### DICOM IMPLEMENTATION CLASS AND VERSION FOR AE STORAGE

|                             |                       |
|-----------------------------|-----------------------|
| Implementation Class UID    | 1.2.410.200001.1.0001 |
| Implementation Version Name | P6A1_0                |

#### 4.2.1.3 Association Initiation Policy

##### 4.2.1.3.1 Activity – Send Images and Structured Reports and Requests Commitment

###### 4.2.1.3.1.1 Description and Sequencing of Activities


A user can select exams or images and request them to be sent to some destination. Each request is forwarded to the job queue and processed individually. When the “Batch” or “Send As You Go” option is active, Stored images and reports will be forwarded to the network job queue for a pre-configured auto-send target destination

automatically. For “Batch” and “Manual” configuration, the system opens an association, sends all images in the study, and closes the association. If “Send As You Go” is selected, the system handles the association with the Storage SCP Server using the following method.

- a. Open an Association when the first image is acquired, and keep association open until the study is closed.
- b. If an error occurs while sending an image to the server because there is no longer an open association (server timed-out), attempt to re-establish the association.
- c. When the study is closed, close the open association after images remained in that study are sent.

Structured Reports are only sent over a separate association at End Exam For.

If the remote AE is configured as an archive device, the Storage AE will, after all images and reports have been sent, transmit Storage Commitment request (N-ACTION) over a separate Association. The Storage AE can only receive an N-EVENT-REPORT request in a subsequent association initiated by the SCP.



**Figure 4.2-1**  
**SEQUENCING OF ACTIVITY - SEND IMAGES**

A possible sequence of interactions between the Storage AE and an Image Manager (e.g. a storage or archive device supporting the Storage and Storage Commitment SOP Classes as an SCP) is illustrated in the figure above.

NOTE: The N-EVENT-REPORT must be sent over a separate association initiated by the Image Manager. (See Section 4.2.1.4)

#### 4.2.1.3.1.2 Proposed Presentation Contexts

SONOACE X8 is capable of proposing the Presentation Contexts shown in the following table.

**Table 4.2-7  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND IMAGES**

| Presentation Context Table              |                                   |  |  |      |              |
|---|-----------------------------------|--|--|------|--------------|
| Abstract Syntax                         |                                   | Transfer Syntax  |  | Role | Ext.<br>Neg. |
| Name                                    | UID                               | Name List  | UID List                                 |      |              |
| Ultrasound Image Storage                | 1.2.840.10008.5.<br>1.4.1.1.6.1   | Implicit VR Little Endian                              | 1.2.840.10008.1.2                        | SCU  | None         |
| Ultrasound Multi-frame Image Storage    | 1.2.840.10008.5.<br>1.4.1.1.3.1   | JPEG Lossy Baseline                                    | 1.2.840.10008.1.2.4<br>.50               | SCU  | None         |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.<br>1.4.1.1.88.33 | Implicit VR Little Endian                              | 1.2.840.10008.1.2                        | SCU  | None         |
| Storage Commitment Push Model           | 1.2.840.10008.1.<br>20.1          | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None         |
| Verification                            | 1.2.840.10008.1.<br>1             | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None         |

Presentation Contexts for Ultrasound Image Storage and Ultrasound Multi-frame Image Storage will be proposed for the “Storage” device configured in Setup/DICOM.

A Presentation Context for Comprehensive Structured Report Storage will be proposed for the “Storage SR” device configured in Setup/DICOM.

A Presentation Context for Storage Commitment Push Model will be proposed for the “SC” device configured in Setup/DICOM.

A Presentation Context for Verification will be proposed when a user press the “Test” button for a configured device.

#### 4.2.1.3.1.3 SOP Specific Conformance Image & Comprehensive Structured Report Storage SOP Classes

All Image & Structured Report Storage SOP Classes supported by the Storage AE exhibit the same behavior, except where stated, and are described together in this section.

**Table 4.2-8  
STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning                   | Error Code             | Behavior   |
|----------------|-----------------------------------|------------------------|--|
| Success        | Success                           | 0000                   | The SCP has successfully stored the SOP Instance. If all SOP Instances succeed, the job is marked as complete. |
| Refused        | Out of Resources                  | A700-A7FF              | The association is aborted using A-ABORT and the send job is marked as failed. The status is logged.           |
| Error          | Data Set does not match SOP Class | A900-A9FF              | Same as "Refused" above.   |
| Error          | Cannot Understand                 | C000-CFFF              | Same as "Refused" above.   |
| Warning        | Coercion of Data Elements         | B000                   | Image transmission is considered successful.   |
| Warning        | Data Set does not match SOP Class | B007                   | Same as "Warning" above.   |
| Warning        | Elements Discards                 | B006                   | Same as "Warning" above.   |
| *              | *                                 | Any other status code. | Same as "Refused" above.   |

The Behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-9  
STORAGE COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. |
| Association aborted by the SCP or network layers | The Send job is marked as failed.  |

A failed send job can be restarted by user interaction. The system can be configured to automatically resend failed

jobs if a transient status code is received. The delay between resending failed jobs and the number of retries is also configurable.

#### 4.2.1.3.1.4 SOP Specific Conformance for Storage Commitment SOP Class

##### 4.2.1.3.1.4.1 Storage Commitment Operations (N-ACTION)

The Storage AE will request storage commitment for the configured device for instances of the Ultrasound Image, Ultrasound Multi-frame Image and Structured Report Storage SOP Classes.

The Storage AE will consider Storage Commitment failed if no N-EVENT-REPORT is received for a Transaction UID within a configurable time period after receiving a successful N-ACTION response (duration of applicability for a Transaction UID).

The Storage AE does not send the optional Storage Media FileSet ID & UID Attributes or the Referenced Study Component Sequence Attribute in the N-ACTION

The Behavior of Storage AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 4.2-10**

**STORAGE COMMITMENT N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code             | Behavior   |
|----------------|-----------------|------------------------|--|
| Success        | Success         | 0000                   | The request for storage comment is considered successfully sent. The system waits for the association of the N-Event-Report. |
| *              | *               | Any other status code. | The Association is aborted using A-Abort and the request for storage comment is marked as failed                             |

The behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-11**

**STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the storage commitment job is marked as failed. |
| Association aborted by the SCP or network layers | The storage commitment job is marked as failed.  |

#### 4.2.1.3.1.4.2 Storage Commitment Notification (N-EVENT-REPORT)

The Storage AE is capable of receiving an N-EVENT-REPORT notification if it has successfully negotiated a Presentation Context for the Storage Commitment Push Model.

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be cancelled.

The behavior of Storage AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below.

**Table 4.2-12  
STORAGE COMMITMENT N-EVENT-REPORT BEHAVIOR**

| Event Type Name                                       | Event Type ID | Behavior  |
|---|---------------|---|
| Storage Commitment Request Successful                 | 1             | The commit status is set to “Y” for each exam in the exam list. Auto deletion for committed exam is not supported.  |
| Storage Commitment Request Complete – Failures Exists | 2             | The commit status is set to “N” for each exam in the exam list. The Referenced SOP Instances under Failed SOP Sequence (0008, 1198) are logged. A send job that failed storage commitment will not be automatically restarted but can be restarted by user interaction. |

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the Table below.

**Table 4.2-13  
STORAGE COMMITMENT N-EVENT-REPORT RESPONSE STATUS REASONS**

| Service Status | Further Meaning        | Error Code | Behavior   |
|----------------|------------------------|------------|--|
| Success        | Success                | 0000       | The Storage commitment result has been successfully received.  |
| Failure        | Unrecognized Operation | 0211H      | The Transaction UID in the N_EVENT_REPORT request is not (was never issued within an N_ACTION request) |
| Failure        | No Such Event Type     | 0113H      | An invalid Event Type ID was supplied in the   |

|         |                    |       |  |
|---------|--------------------|-------|--|
|         |                    |       | N_EVENT_REPORT request   |
| Failure | Processing Failure | 0110H | An internal error occurred during processing of the N_EVENT_REPORT |

#### 4.2.1.3.1.5 SOP Specific Conformance for Verification

The Behavior when encountering status codes in a C-ECHO response is summarized in the Table below:

**Table 4.2-14  
VERIFICATION C-ECHO RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code            | Behavior                               |
|----------------|-----------------|-----------------------|--|
| Success        | Success         | 0000                  | Verification Status is set to 'Normal' |
| *              | *               | Any other status code | Verification Status is set to 'Failed' |

The Behavior of Storage AE during communication failure is summarized in the Table below:

**Table 4.2-15  
VERIFICATION COMMUNICATION FAILURE BEHAVIOR**


| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the verification job is marked as failed. |
| Association aborted by the SCP or network layers | The verification job is marked as failed.  |

#### 4.2.1.4 Association Acceptance Policy

##### 4.2.1.4.1 Activity – Receive Storage Commitment Response

##### 4.2.1.4.1.1 Description and Sequence of Activities

The Storage AE will accept associations in order to receive responses to a Storage Commitment Request.



**Figure 4.2-2**  
**SEQUENCING OF ACTIVITY - RECEIVE STORAGE COMMITMENT RESPONSE**

A possible sequence of interactions between the Storage AE and an Image Manager (e.g. a storage or archive device supporting Storage Commitment SOP Classes as an SCP) is illustrated in the Figure above:

1. The Image Manager opens a new association with the Storage AE.
2. The Image Manager sends an N-EVENT-REPORT request notifying the Storage AE of the status of a previous Storage Commitment Request. The Storage AE replies with an N-EVENT-REPORT response confirming receipt.
3. The Image Manager closes the association with the Storage AE.

**4.2.1.4.1.2 Accepted Presentation Contexts**

The Storage AE will accept Presentation Contexts as shown in the Table below.

**Table 4.2-16**  
**ACCEPTABLE PRESENTATION CONTEXTS FOR ACTIVITY**  
**RECEIVE STORAGE COMMITMENT RESPONSE**

| Presentation Context Table |     |                 |          |      |           |
|----------------------------|-----|-----------------|----------|------|-----------|
| Abstract Syntax            |     | Transfer Syntax |          | Role | Ext. Neg. |
| Name                       | UID | Name List       | UID List |      |           |
|                            |     |                 |          |      |           |

|                               |                      |  |  |     |      |
|-------------------------------|----------------------|--|--|-----|------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU | None |
| Verification                  | 1.2.840.10008.1.1    | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP | None |

#### 4.2.1.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

##### 4.2.1.4.1.3.1 Storage Commitment Notifications (N-EVENT-REPORT)

Upon receipt of an N-EVENT-REPORT the timer associated with the Transaction UID will be cancelled.

The behavior of Storage AE when receiving Event Types within the N-EVENT\_REPORT is summarized in Table 4.2-12.

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in Table 4.2-13.

##### 4.2.1.4.1.4 SOP Specific Conformance for Verification SOP Class

The Storage AE provides standard conformance to the Verification SOP Class as an SCP. If the C-ECHO request was successfully received, a 0000 (Success) status code will be returned in the C-ECHO response.

## 4.2.2 Workflow Application Entity Specification

### 4.2.2.1 SOP Classes

SONOACE X8 provides Standard Conformance to the following SOP Classes:

**Table 4.2-17**  
**SOP CLASSES FOR AE WORKFLOW**

| SOP Classes                                | SOP Class UID           | SCU | SCP |
|--|-------------------------|-----|-----|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31  | Yes | No  |
| Modality Performed Procedure Step          | 1.2.840.10008.3.1.2.3.3 | Yes | No  |

## 4.2.2.2 Association Establishment Policy

### 4.2.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed.

**Table 4.2-18**

**DICOM APPLICATION CONTEXT FOR AE WORKFLOW**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

### 4.2.2.2.2 Number of Associations

SONOACE X8 initiates one Association at a time for a Worklist request.

**Table 4.2-19**

**NUMBER OF ASSOCIATIONS INITIATED FOR AE WORKFLOW**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

### 4.2.2.2.3 Asynchronous Nature

SONOACE X8 does not support asynchronous communications (multiple outstanding transactions over a single Association)

**Table 4.2-20**

**ASYNCHRONOUS NATURE AS A SCU FOR AE WORKFLOW**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

### 4.2.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-21**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE WORKFLOW**

|                             |                       |
|-----------------------------|-----------------------|
| Implementation Class UID    | 1.2.410.200001.1.0001 |
| Implementation Version Name | P6A1_0                |

### **4.2.2.3 Association Initiation Policy**

#### **4.2.2.3.1 Activity – Worklist Update**

##### **4.2.2.3.1.1 Description and Sequencing of Activities**


The request for a Worklist Update is initiated by user interaction or automatically at specific time intervals, configurable by the user.

The interactive Worklist Query will display a dialog for entering data as search criteria. When the Query is started on your request, only the data from the dialog will be inserted as matching keys into the query.

With automated worklist queries the SONOACE X8 always requests all items for a Scheduled Procedure Step Start Date (actual date), Modality (US) and Scheduled Station AE Title.

Upon initiation of the request, the SONOACE X8 will build an Identifier for the C-FIND request, will initiate an Association to send the request and will wait for Worklist responses. After retrieval of all responses, SONOACE X8 will access the local database to add patient demographic data. The results will be displayed in a separate list, which will be cleared with the next worklist update.

SONOACE X8 will initiate an Association in order to issue a C-FIND request according to the Modality Worklist Information Model.



**Figure 4.2-3**  
**SEQUENCING OF ACTIVITY - WORKLIST UPDATE**

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as an SCP) is illustrated in the figure above:

**4.2.2.3.1.2 Proposed Presentation Contexts**

SONOACE X8 will propose Presentation Contexts as shown in the following table:

**Table 4.2-22**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY WORKLIST UPDATE**

| Presentation Context Table                 |                |                           |                      |      |           |
|--|----------------|---------------------------|----------------------|------|-----------|
| Abstract Syntax                            |                | Transfer Syntax           |                      | Role | Ext. Neg. |
| Name                                       | UID            | Name List                 | UID List             |      |           |
| Modality Worklist Information Model - FIND | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2    | SCU  | None      |
|  | 5.1.4.31       | Explicit VR Little Endian | 1.2.840.10008. 1.2.1 |      |           |

#### 4.2.2.3.1.3 SOP Specific Conformance for Modality Worklist

The behavior of SONOACE X8 when encountering status codes in a Modality Worklist C-FIND response is summarized in the Table below. If any other SCP response status than “Success” or “Pending” is received by SONOACE X8, a message “Query failed” will appear on the user interface.

**Table 4.2-23  
MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Error Code</b>            | <b>Behavior</b>  |
|-----------------------|---|------------------------------|--|
| Success               | Matching is complete  | 0000                         | The SCP has Completed the operation successfully.                                |
| Pending               | Matches are continuing  | FF00                         | Continue.  |
| Pending               | Matches are continuing -<br>Warning that one or more<br>Optional Keys were not<br>supported | FF01                         | Continue.  |
| *                     | *   | Any other<br>status<br>code. | The Association is aborted using A-Abort and the<br>Worklist is marked as failed |

The behavior of SONOACE X8 during communication failure is summarized in the Table below.

**Table 4.2-24  
MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>   |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and<br>the worklist query is marked as failed. |
| Association aborted by the SCP or network layers | The Worklist query is marked as failed.   |

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The Table below provides a description of the SONOACE X8 Worklist Request Identifier and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

Requested return attributes not supported by the SCP are set to have no value. Non-matching responses returned by the SCP due to unsupported optional matching keys are ignored. No attempt is made to filter out possible duplicate entries.

**Table 4.2-25  
WORKLIST REQUEST IDENTIFIER**

| Module Name<br>Attribute Name           | Tag       | VR | M   | R | Q | D | IOD |
|---|-----------|----|-----|---|---|---|-----|
| Scheduled Procedure Step                |           |    |     |   |   |   |     |
| Scheduled Procedure Step Sequence       | 0040,0100 | SQ |     | x |   |   |     |
| > Scheduled Station AET                 | 0040,0001 | AE | (S) | x | x |   |     |
| > Scheduled Procedure Step Start Date   | 0040,0002 | DA | S,R | x | x | x |     |
| > Scheduled Procedure Step Start Time   | 0040,0003 | TM |     | x |   | x |     |
| > Modality                              | 0008,0060 | CS | S   | x | x |   |     |
| > Scheduled Performing Physician's Name | 0040,0006 | PN |     | x |   |   |     |
| > Scheduled Procedure Step Description  | 0040,0007 | LO |     | x |   | x | x   |
| > Scheduled Station Name                | 0040,0010 | SH |     | x |   |   |     |
| > Scheduled Procedure Step Location     | 0040,0011 | SH |     | x |   |   |     |
| > Scheduled Protocol Code Sequence      | 0040,0008 | SQ |     | x |   |   | x   |
| > Scheduled Procedure Step ID           | 0040,0009 | SH |     | x |   |   | x   |
| Requested Procedure                     |           |    |     |   |   |   |     |
| Requested Procedure ID                  | 0040,1001 | SH |     | x |   | x | x   |
| Requested Procedure Description         | 0032,1060 | LO |     | x |   |   |     |
| Study Instance UID                      | 0020,000D | UI |     | x |   |   | x   |
| Referenced Study Sequence               | 0008,1110 | SQ |     | x |   |   | x   |
| Requested Procedure Code Sequence       | 0032,1064 | SQ |     | x |   |   | x   |
| Imaging Service Request                 |           |    |     |   |   |   |     |
| Accession Number                        | 0008,0050 | SH |     | x |   | x | x   |
| Requesting Physician                    | 0032,1032 | PN |     | x |   |   |     |
| Referring Physician's Name              | 0008,0090 | PN |     | x |   |   | x   |
| Visit Status                            |           |    |     |   |   |   |     |
| Current Patient Location                | 0038,0300 | LO |     | x |   |   |     |
| Patient Identification                  |           |    |     |   |   |   |     |
| Patient's Name                          | 0010,0010 | PN |     | x |   | x | x   |
| Patient ID                              | 0010,0020 | LO |     | x |   | x | x   |
| Patient Demographic                     |           |    |     |   |   |   |     |

|                      |           |    |   |   |   |
|----------------------|-----------|----|---|---|---|
| Patient's Birth Date | 0010,0030 | DA | x | x | x |
| Patient's Sex        | 0010,0040 | CS | x | x | x |
| Patient's Size       | 0010,1020 | DS | x | x | x |
| Patient's Weight     | 0010,1030 | DS | x | x | x |

The above table should read as follows:

- Module Name: The Name of the associated module for supported worklist attributes.
- Attribute Name: Attributes supported to build an SONOACE X8 Worklist Request Identifier.
- Tag: DICOM tag for this attribute.
- VR: DICOM VR for this attribute.
- M: Matching keys for (automatic) Worklist Update. An "S" indicates that SONOACE X8 supplies an attribute value for Single Value Matching or additional specific tags indicated by "(S)"; an "R" will indicate Range Matching.
- R: Return keys. An "X" will indicate that SONOACE X8 will supply this attribute as Return Key with zero length for Universal Matching.
- Q: Interactive Query Key. An "X" will indicate that SONOACE X8 will supply this attribute as matching key, if entered in the Setup Dialog.
- D: Displayed keys. An "X" indicates that this worklist attribute is displayed to the user during a patient registration dialog.
- IOD: An "X" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

#### 4.2.2.3.2 Activity – Acquire Images

##### 4.2.2.3.2.1 Description and Sequencing of Activities

An Association to the configured MPPS SCP system is established immediately after the first image is acquired to send the MPPS N-Create message.


The "End Exam" button causes a "COMPLETED" message. An exam for which an MPPS instance is sent with a state of "COMPLETED" can no longer be updated.

The SONOACE X8 will support creation of "unscheduled cases" by allowing MPPS Instances to be communicated for locally registered Patients.

The SONOACE X8 only supports a 1-to-1 relationship between Scheduled and Performed Procedure Steps.

SONOACE X8 will initiate an Association to issue an:

- N-CREATE request according to the CREATE Modality Performed Procedure Step SOP Instance operation, or an:
- N-SET request to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.



**Figure 4.2-4**  
**SEQUENCING OF ACTIVITY - ACQUIRE IMAGES**

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the MPPS SOP Class as an SCP) is illustrated in the figure above:

#### 4.2.2.3.2 Proposed Presentation Contexts

SONOACE X8 will propose Presentation Contexts as shown in the following table:

**Table 4.2-26  
PROPOSED PRESENTATION CONTEXTS FOR REAL-WORLD ACTIVITY ACQUIRE IMAGES**

| Presentation Context Table |                |                           |                      |      |           |
|----------------------------|----------------|---------------------------|----------------------|------|-----------|
| Abstract Syntax            |                | Transfer Syntax           |                      | Role | Ext. Neg. |
| Name                       | UID            | Name List                 | UID List             |      |           |
| Modality Performed         | 1.2.840.10008. | Implicit VR Little Endian | 1.2.840.10008.1.2    | SCU  | None      |
| Procedure Step             | 3.1.2.3.3      | Explicit VR Little Endian | 1.2.840.10008. 1.2.1 |      |           |

#### 4.2.2.3.2.3 SOP Specific Conformance for MPPS

The behavior of SONOACE X8 when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in the Table below. If any other SCP response status than “Success” or “Warning” is received by SONOACE X8, a message “MPPS failed” will appear on the user interface.

**Table 4.2-27  
MPPS N-CREATE / N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning              | Error Code             | Behavior  |
|----------------|------------------------------|------------------------|---|
| Success        | Success                      | 0000                   | The SCP has Completed the operation successfully.                         |
| Warning        | Attribute Value Out of Range | 0116H                  | The MPPS Operation is considered successful.                              |
| *              | *                            | Any other status code. | The Association is aborted using A-Abort and the MPPS is marked as failed |

The behavior of SONOACE X8 during communication failure is summarized in the table below:

**Table 4.2-28**

### MPPS COMMUNICATION FAILURE BEHAVIOR

| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the MPPS job is marked as failed. |
| Association aborted by the SCP or network layers | The MPPS job is marked as failed.  |

Table 4.2-29 provides a description of the MPPS N-CREATE and N-SET request identifiers send by SONOACE X8. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent.

**Table 4.2-29**  
**MPPS N-CREATE / N-SET REQUEST IDENTIFIER**

| Attribute Name                               | Tag       | VR | N-CREATE                        | N-SET |
|--|-----------|----|---------------------------------|-------|
| Specific Character Set                       | 0008,0005 | CS | "ISO_IR 100"                    |       |
| <b>Performed Procedure Step Relationship</b> |           |    |                                 |       |
| Scheduled Step Attribute Sequence            | 0040,0270 | SQ |                                 |       |
| > Study Instance UID                         | 0020,000D | UI | From MWL or generated by device |       |
| > Referenced Study Sequence                  | 0008,1110 | SQ | From MWL                        |       |
| >> Referenced SOP Class UID                  | 0008,1150 | UI | From MWL                        |       |
| >> Referenced SOP Instance UID               | 0008,1155 | UI | From MWL                        |       |
| > Accession Number                           | 0008,0050 | SH | From MWL or user input          |       |
| > Requested Procedure ID                     | 0040,1001 | SH | From MWL                        |       |
| > Requested Procedure Description            | 0032,1060 | LO | From MWL                        |       |
| > Scheduled Procedure Step ID                | 0040,0009 | SH | From MWL                        |       |
| > Scheduled Procedure Step Description       | 0040,0007 | LO | From MWL                        |       |
| > Scheduled Protocol Code Sequence           | 0040,0008 | SQ | From MWL                        |       |
| >> Code Value                                | 0008,0100 | SH | From MWL                        |       |
| >> Coding Scheme Designator                  | 0008,0102 | SH | From MWL                        |       |
| >> Coding Scheme Version                     | 0008,0103 | SH | From MWL                        |       |

|   |           |    |  |  |
|---|-----------|----|--|--|
| >> Code Meaning                             | 0008,0104 | LO | From MWL   |  |
| Patient's Name                              | 0010,0010 | PN | From MWL or user input                                   |  |
| Patient ID                                  | 0010,0020 | LO | From MWL or user input                                   |  |
| Patient's Birth Date                        | 0010,0030 | DA | From MWL or user input                                   |  |
| Patient's Sex                               | 0010,0040 | CS | From MWL or user input                                   |  |
| Referenced Patient Sequence                 | 0008,1120 | SQ | Zero length  |  |
| > Referenced SOP Class UID                  | 0008,1150 | UI | Zero length  |  |
| > Referenced Instance UID                   | 0008,1155 | UI | Zero length  |  |
| <b>Performed Procedure Step Information</b> |           |    |  |  |
| Performed Procedure Step ID                 | 0040,0253 | SH | Generated by device<br>(Study Date + Study<br>Time)      |  |
| Performed Station AE Title                  | 0040,0241 | AE | From Modality Setup                                      |  |
| Performed Station Name                      | 0040,0242 | SH | Zero length  |  |
| Performed Location                          | 0040,0243 | SH | Zero length  |  |
| Performed Procedure Step<br>Start Date      | 0040,0244 | DA | Actual Start Date  |  |
| Performed Procedure Step<br>Start Time      | 0040,0245 | TM | Actual Start Time  |  |
| Performed Procedure Step<br>Status          | 0040,0252 | CS | "IN PROGRESS"  | "DISCONTINUED" or<br>"COMPLETED"                         |
| Performed Procedure Step<br>Description     | 0040,0254 | LO | From MWL or user input<br>(Same as Study<br>Description) | From MWL or user input<br>(Same as Study<br>Description) |
| Performed Procedure Type<br>Description     | 0040,0255 | LO | Zero length  | Zero length  |
| Procedure Code Sequence                     | 0008,1032 | SQ | From MWL   | From MWL   |
| > Code Value                                | 0008,0100 | SH | From MWL   | From MWL   |
| > Coding Scheme Designator                  | 0008,0102 | SH | From MWL   | From MWL   |
| > Coding Scheme Version                     | 0008,0103 | SH | From MWL   | From MWL   |
| > Code Meaning                              | 0008,0104 | LO | From MWL   | From MWL   |
| Performed Procedure Step End<br>Date        | 0040,0250 | DA | Zero length  | Actual End Date  |

|  |           |    |  |   |
|--|-----------|----|--|---|
| Performed Procedure Step End Time                      | 0040,0251 | TM | Zero length                                      | Actual End Time   |
| <b>Image Acquisition Results</b>                       |           |    |  |   |
| Modality   | 0008,0060 | CS | "US"   |   |
| Study ID   | 0020,0010 | SH | generated by device<br>(Study Date + Study Time) |   |
| Performed Protocol Code Sequence                       | 0040,0260 | SQ | Zero length                                      |   |
| Performed Series Sequence                              | 0040,0340 | SQ | Zero length                                      | One or more items   |
| > Performed Physician's Name                           | 0008,1050 | PN |  | From MWL  |
| > Protocol Name  | 0018,1030 | LO |  | "FreeForm" or staged protocol name (In case of Stress echo) |
| > Operator's Name                                      | 0008,1070 | PN |  | From user input   |
| > Series Instance UID                                  | 0020,000E | UI |  | generated by device   |
| > Series Description                                   | 0008,103E | LO |  | Zero length   |
| > Retrieve AE Title                                    | 0008,0054 | AE |  | Zero length   |
| > Referenced Image Sequence                            | 0008,1140 | SQ |  | From Modality   |
| >> Referenced SOP Class UID                            | 0008,1150 | UI |  | From Modality   |
| >> Referenced SOP Instance UID                         | 0008,1155 | UI |  | From Modality   |
| > Referenced Non-Image Composite SOP Instance Sequence | 0040,0220 | SQ |  | From Modality   |
| >> Referenced SOP Class UID                            | 0008,1150 | UI |  | From Modality   |
| >> Referenced SOP Instance UID                         | 0008,1155 | UI |  | From Modality   |

#### 4.2.2.4 Association Acceptance Policy

The Workflow Application Entity does not accept Associations.

## 4.2.3 Hardcopy Application Entity Specification

### 4.2.3.1 SOP Classes

SONOACE X8 provides Standard Conformance to the following SOP Classes:

**Table 4.2-30**  
**SOP CLASSES FOR AE HARDCOPY**

| SOP Classes                           | SOP Class UID          | SCU | SCP |
|---------------------------------------|------------------------|-----|-----|
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9  | Yes | No  |
| Basic Color Print Management Meta     | 1.2.840.10008.5.1.1.18 | Yes | No  |

### 4.2.3.2 Association Policies

#### 4.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-31**  
**DICOM APPLICATION CONTEXT FOR AE HARDCOPY**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.3.2.2 Number of Association

SONOACE X8 can initiate one or more Associations at a time for each destination to which a transfer request is being processed in the active job queue list.

**Table 4.2-32**  
**NUMBER OF ASSOCIATIONS INITIATED FOR AE HARDCOPY**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | Unlimited (number of configured hardcopy devices) |
|---|---|

#### 4.2.3.2.3 Asynchronous Nature

SONOACE X8 does not support asynchronous communications (multiple outstanding transactions over a single Association)

**Table 4.2-33**  
**ASYNCHRONOUS NATURE AS A SCU FOR AE HARDCOPY**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

#### 4.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-34**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR AE HARDCOPY**


|                             |                       |
|-----------------------------|-----------------------|
| Implementation Class UID    | 1.2.410.200001.1.0001 |
| Implementation Version Name | P6A1_0                |

#### 4.2.3.3 Association Initiation Policy

##### 4.2.3.3.1 Activity – Film Images

##### 4.2.3.3.1.1 Description and Sequencing of Activities

A user composes images onto film sheets and requests them to be sent to a specific hardcopy device. The user can select the desired film format and number of copies. Each print-job is forwarded to the job queue and processed individually.



**Figure 4.2-5**

## SEQUENCING OF ACTIVITY - FILM IMAGES

A typical sequence of DIMSE messages sent over an association between Hardcopy AE and a Printer is illustrated in the Figure above:

Association Initiation Policies for “Batch”, “Send As You Go” and “Manual” Mode are equal to the Sending images’ of the Storage Application Entity. (See 4.2.1.3.1.1)

Status of the print-job is reported through the job control interface. One or more job can be active at a time for each separate hardcopy device. If any response from the remote Application contains a status other than Success or Warning, the Association is aborted and the related job is switched to a failed state. It can be restarted any time by user interaction or, if configured, by automated retry.

### 4.2.3.3.1.2 Proposed Presentation Contexts

SONOACE X8 is capable of proposing the Presentation Contexts shown in the Table below:

**Table 4.2-35  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY FILM IMAGES**

| Presentation Context Table            |                        |  |  |      |           |
|---------------------------------------|------------------------|--|--|------|-----------|
| Abstract Syntax                       |                        | Transfer Syntax  |  | Role | Ext. Neg. |
| Name                                  | UID                    | Name List  | UID List                                 |      |           |
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9  | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |
| Basic Color Print Management Meta     | 1.2.840.10008.5.1.1.18 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |

### 4.2.3.3.1.3 Common SOP Specific Conformance for all Print SOP Classes

The general behavior of Hardcopy AE during communication failure is summarized in the table below. This behavior is common for all SOP Classes supported by Hardcopy AE.

**Table 4.2-36  
HARDCOPY COMMUNICATION FAILURE BEHAVIOR**

| Exception | Behavior  |
|-----------|---|
| Timeout   | The Association is aborted using A-ABORT and the print job is marked as failed. |

|  |                                    |
|--|------------------------------------|
| Association aborted by the SCP or network layers | The print job is marked as failed. |
|--|------------------------------------|

#### 4.2.3.3.1.4 SOP Specific Conformance for the Film Session SOP Class

Hardcopy AE supports the following DIMSE operations for the Film Session SOP Class:

- N-CREATE
- N-DELETE

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.3.3.1.4.1 Film Session SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below:

**Table 4.2-37**  
**FILM SESSION SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| Attribute Name   | Tag       | VR | Value   | Presence of Value | Source |
|------------------|-----------|----|---|-------------------|--------|
| Number of Copies | 2000,0010 | IS | 1..99   | ALWAYS            | USER   |
| Print Priority   | 2000,0020 | CS | HIGH, MED or LOW  | ALWAYS            | USER   |
| Medium Type      | 2000,0030 | CS | PAPER, CLEAR FILM, BLUE FILM, MAMMO CLEAR FILM or MAMMO BLUE FILM | ALWAYS            | USER   |
| Film Destination | 2000,0040 | CS | MAGAZINE or PROCESSOR   | ALWAYS            | USER   |

The Behavior of Hardcopy AE when encountering status codes in an N-CREATE response is summarized in the table below:

**Table 4.2-38**  
**FILM SESSION SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning     | Error Code | Behavior  |
|----------------|---------------------|------------|---|
| Success        | Success             | 0000       | The SCP has Completed the operation successfully. |
| Warning        | Attribute Value Out | 0116H      | System continues operations.                      |

|         |                      |                        |  |
|---------|----------------------|------------------------|--|
|         | of Range             |                        |  |
| Warning | Attribute List Error | 0107H                  | Same as above  |
| *       | *                    | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.4.2 Film Session SOP Class Operations (N-DELETE)

The behavior of Hardcopy AE when encountering status codes in an N-DELETE response is summarized in the Table below:

**Table 4.2-39**

#### **PRINTER SOP CLASS N-DELETE RESONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Error Code</b>      |  |
|-----------------------|------------------------|------------------------|--|
| Success               | Success                | 0000                   | The SCP has Completed the operation successfully.                              |
| *                     | *                      | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.5 SOP Specific Conformance for the Film Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

#### 4.2.3.3.1.5.1 Film Box SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the table below:

**Table 4.2-40**

**FILM BOX SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| <b>Attribute Name</b>            | <b>Tag</b> | <b>VR</b> | <b>Value</b>  | <b>Presence of Value</b> | <b>Source</b> |
|----------------------------------|------------|-----------|---|--------------------------|---------------|
| Image Display Format             | 2010,0010  | ST        | "STANDARD\1, 1" ,<br>"STANDARD\1, 2" ,<br>"STANDARD\2, 2" ,<br>"STANDARD\2, 3" ,<br>"STANDARD\3, 3" ,<br>"STANDARD\3, 4" ,<br>"STANDARD\3, 5" ,<br>"STANDARD\4, 4" ,<br>"STANDARD\4, 5" or<br>"STANDARD\4, 6" | ALWAYS                   | USER          |
| Referenced Film Session Sequence | 2010.0500  | SQ        |   | ALWAYS                   | AUTO          |
| > Referenced SOP Class UID       | 0008,1150  | UI        | 1.2.840.10008.5.1.1.1   | ALWAYS                   | AUTO          |
| > Referenced SOP Instance UID    | 0008,1155  | UI        | From created Film Session SOP Instance  | ALWAYS                   | AUTO          |
| Film Orientation                 | 2010,0040  | CS        | PORTRAIT or LANDSCAPE   | ALWAYS                   | USER          |
| Film Size ID                     | 2010,0050  | CS        | 8INX10IN, 8_5INX11IN,<br>10INX12IN, 10INX14IN,<br>11INX14IN, 11INX17IN,<br>14INX14IN, 14INX17IN,<br>24CMX24CM, 24CMX30CM,<br>A4, A3   | ALWAYS                   | USER          |
| Magnification Type               | 2010,0060  | CS        | REPLICATE, BILINEAR,<br>CUBIC, NONE   | ALWAYS                   | USER          |
| Max Density                      | 2010,0130  | US        | 0 ~   | ANAP                     | USER          |
| Configuration Information        | 2010,0150  | ST        | Values are defined in Print Conformance Statement   | ANAP                     | USER          |
| Smoothing Type                   | 2010,0080  | CS        | Values are defined in Print Conformance Statement   | ANAP                     | USER          |

|                     |           |    |                |        |      |
|---------------------|-----------|----|----------------|--------|------|
| Border Density      | 2010,0100 | CS | BLACK or WHITE | ALWAYS | USER |
| Empty Image Density | 2010,0110 | CS | BLACK or WHITE | ALWAYS | USER |
| Min Density         | 2010,0120 | US | 0 ~            | ANAP   | USER |

The behavior of Hardcopy AE when encountering status codes in an N-CREATE responses is summarized in the table below:

**Table 4.2-41  
FILM BOX SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Error Code</b>      | <b>Behavior</b>  |
|-----------------------|---|------------------------|--|
| Success               | Success   | 0000                   | The SCP has Completed the operation successfully.                              |
| Warning               | Attribute Value Out of Range  | 0116H                  | System continues operations.   |
| Warning               | Attribute List Error  | 0107H                  | Same as above  |
| Warning               | Requested Min Density or Max Density outside of printer's operating range | B605H                  | Same as above  |
| *                     | *   | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.5.2 Film Box SOP Class Operations (N-ACTION)

An N-ACTION Request is issued to instruct the Print SCP to print the contents of the Film Box.

The behavior of Hardcopy AE when encountering status codes in an N-ACTION responses is summarized in the table below:

**Table 4.2-42  
FILM BOX CLASS N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code             | Behavior   |
|----------------|-----------------|------------------------|--|
| Success        | Success         | 0000                   | The SCP has Completed the operation successfully.                              |
| *              | *               | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |

#### 4.2.3.3.1.6 SOP Specific Conformance for the Film Box SOP Class

Hardcopy AE supports the following DIMSE operations for the Image Box SOP Class:

- N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.3.3.1.6.1 Image Box SOP Class Operations (N-SET)

The attributes supplied in an N-SET Request are listed in the Table below:

**Table 4.2-43  
BASIC GRAYSCALE IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES**

| Attribute Name                 | Tag       | VR | Value                                 | Presence of Value | Source |
|--------------------------------|-----------|----|---------------------------------------|-------------------|--------|
| Image Position                 | 2020,0010 | US | 1 .. N (N = Row * Column of Film Box) | ALWAYS            | AUTO   |
| Basic Grayscale Image Sequence | 2020,0110 | SQ |                                       | ALWAYS            | AUTO   |
| > Samples Per Pixel            | 0028,0002 | US | 1                                     | ALWAYS            | AUTO   |
| > Photometric Interpretation   | 0028,0004 | CS | MONOCHROME2                           | ALWAYS            | AUTO   |
| > Rows                         | 0028,0010 | US | Number of Row Pixels of Image         | ALWAYS            | AUTO   |
| > Columns                      | 0028,0011 | US | Number of Column Pixels of Image      | ALWAYS            | AUTO   |
| > Bits Allocated               | 0028,0100 | US | 8                                     | ALWAYS            | AUTO   |
| > Bits Stored                  | 0028,0101 | US | 8                                     | ALWAYS            | AUTO   |

|                        |           |    |                 |        |      |
|------------------------|-----------|----|-----------------|--------|------|
| > High Bit             | 0028,0102 | US | 7               | ALWAYS | AUTO |
| > Pixel Representation | 0028,0103 | US | 0               | ALWAYS | AUTO |
| > Pixel Data           | 7FE0,0010 | OB | Pixels of Image | ALWAYS | AUTO |

**Table 4.2-44**

**BASIC COLOR IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES**

| Attribute Name               | Tag       | VR | Value                                 | Presence of Value | Source |
|------------------------------|-----------|----|---------------------------------------|-------------------|--------|
| Image Position               | 2020,0010 | US | 1 .. N (N = Row * Column of Film Box) | ALWAYS            | AUTO   |
| > Samples Per Pixel          | 0028,0002 | US | 3                                     | ALWAYS            | AUTO   |
| > Photometric Interpretation | 0028,0004 | CS | RGB                                   | ALWAYS            | AUTO   |
| > Planar Configuration       | 0028,0006 | US | 1                                     | ALWAYS            | AUTO   |
| > Rows                       | 0028,0010 | US | Number of Row Pixels of Image         | ALWAYS            | AUTO   |
| > Columns                    | 0028,0011 | US | Number of Column Pixels of Image      | ALWAYS            | AUTO   |
| > Bits Allocated             | 0028,0100 | US | 8                                     | ALWAYS            | AUTO   |
| > Bits Stored                | 0028,0101 | US | 8                                     | ALWAYS            | AUTO   |
| > High Bit                   | 0028,0102 | US | 7                                     | ALWAYS            | AUTO   |
| > Pixel Representation       | 0028,0103 | US | 0                                     | ALWAYS            | AUTO   |
| > Pixel Data                 | 7FE0,0010 | OB | Pixels of Image                       | ALWAYS            | AUTO   |

The behavior of Hardcopy AE when encountering status codes in an N-SET response is summarized in the table below:

**Table 4.2-45**

**IMAGE BOX SOP CLASS N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Error Code | Behavior  |
|----------------|-----------------|------------|---|
| Success        | Success         | 0000       | The SCP has Completed the operation successfully. |

|   |   |                        |  |
|---|---|------------------------|--|
| * | * | Any other status code. | The Association is aborted using A-Abort and the print-job is marked as failed |
|---|---|------------------------|--|

#### 4.2.3.4 Association Acceptance Policy

The Hardcopy Application Entity does not accept Associations.

## 4.3 NETWORK INTERFACE

### 4.3.1 Physical Network Interface

SONOACE X8 supports a single network interface. One of the following physical network interfaces will be available depending on hardware options installed:

**Table 4.3-1**  
**SUPPORTED PHYSICAL NETWORK INTERFACES**

|                   |
|-------------------|
| Ethernet 100baseT |
| Ethernet 10baseT  |

## 4.4 CONFIGURATION

### 4.4.1 AE Title/Presentation Address Mapping

#### 4.4.1.1 Local AE Titles

All local applications use the AE Titles and TCP/IP Ports configured via the Setup/DICOM Menu. All local DICOM services use the same AE Title. The system listens for Verification requests and Commitment reports on the configured Port.

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

The AE Title, host names and port numbers of remote applications are configured using the SONOACE X8 Setup/DICOM Menu.

#### 4.4.1.2.1 Storage

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Image Storage SCPs. Multiple remote Image Storage SCPs can be defined.

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Structured Report Storage SCP. Only a single remote Structured Report

Storage SCP can be defined.

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Storage Commitment SCP. Only a single remote Storage Commitment SCP can be defined and only one Image Storage SCP can be assigned for Storage Commitment.

#### 4.4.1.2.2 Workflow

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Modality Worklist SCP. Only a single remote Modality Worklist SCP can be defined.

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote MPPS SCP. Only a single remote MPPS SCP can be defined.

#### 4.4.1.2.3 Hardcopy

The Add button on the SONOACE X8 Setup/DICOM Menu must be used to set the AE Titles, port-numbers, IP addresses and capabilities for the remote Print SCPs. Multiple remote Print SCPs can be defined.

### 4.4.2 Parameters

A number of parameters related to acquisition and general operation can be configured using the Setup/DICOM Menu. The Table below only shows those configuration parameters relevant to DICOM communications. See the SONOACE X8 Manual for details on general configuration capabilities.

**Table 4.4-1  
CONFIGURATION PARAMETERS TABLE**


| <b>Parameter</b>                 | <b>Configurable<br/>(Yes/No)</b> | <b>Default Value</b> |
|----------------------------------|----------------------------------|----------------------|
| Local System Parameters          |                                  |                      |
| AE Title (Local System AE Title) | Yes                              | “Set AE Title”       |
| Station Name                     | Yes                              | “Set Station Name”   |
| Port No. (Local Port Number)     | Yes                              | 104                  |

| <b>Service Common Parameters</b>                       |     |              |
|--|-----|--------------|
| Retry Interval   | Yes | 30 Sec.      |
| Connect Timeout  | Yes | 15 Sec.      |
| Maximum Retires  | Yes | 1            |
| <b>Storage Parameters</b>                              |     |              |
| Transfer Mode  | Yes | "Batch"      |
| Include 3D Volume                                      | Yes | Checked      |
| Window Center (VOI LUT)                                | Yes | 128          |
| Window Width (VOI LUT)                                 | Yes | 256          |
| <b>Storage Commitment Parameters</b>                   |     |              |
| Associated Storage Server                              | Yes | None         |
| <b>Worklist Modality Parameters</b>                    |     |              |
| Delay between automatic Worklist Updates               | Yes | 5 Min.       |
| Query Worklist for specific Scheduled Station AE Title | Yes | Any          |
| Query Worklist for specific Scheduled Modality Value   | No. | "US" fixed   |
| Query Worklist for Specific Start Date                 | Yes | Today        |
| <b>Print Parameters</b>                                |     |              |
| Transfer Mode  | Yes | "Batch"      |
| Color  | Yes | "Grayscale"  |
| Medium Type  | Yes | "PAPER"      |
| Format   | Yes | 1x1          |
| Film Size  | Yes | 8 IN X 10 IN |
| Orientation  | Yes | "PORTRAIT"   |
| Destination  | Yes | "MAGAZINE"   |
| Magnification  | Yes | "REPLICATE"  |
| Smoothing Type   | Yes | Blank        |
| Border Density   | Yes | "BLACK"      |
| Empty Density  | Yes | "BLACK"      |
| Priority   | Yes | "HIGH"       |
| Min Density  | Yes | Blank        |
| Max Density  | Yes | Blank        |
| Copies   | Yes | 1            |
| Configuration Info                                     | Yes | Blank        |

## 5 MEDIA INTERCHANGE

### 5.1 IMPLEMENTATION MODEL

#### 5.1.1 Application Data Flow



**Figure 5.1-1**

#### **APPLICATION DATA FLOW DIAGRAM FOR MEDIA STORAGE**

- The Offline-Media Application Entity exports images and Structured Report to a DVD Storage medium. It is associated with the local real-world activity "Export to DVD RW", "Export to DVD RW" is performed upon user request for selected studies.

#### 5.1.2 Functional Definition of AEs

##### 5.1.2.1 Functional Definition of Offline-Media Application Entity

Activation of the "Export to DVD RW" menu entry will pass the currently selected studies to the Offline-Media Application Entity. The SOP Instances associated with the selection will be collected into one or more export jobs. The contents of each export job will be written to a single DVD media.

#### 5.1.3 Sequencing of Real-World Activities

At least one study must exist and be selected before the Offline-Media Application Entity can be invoked. The operator can insert a new DVD media at any time before or after invocation of the Offline-Media Application Entity. If no DVD media is available the export job can be cancelled immediately.

#### 5.1.4 File Meta Information Options

The implementation written to the File Meta Header in each file is:

**Table 5.1-1**

#### **DICOM IMPLEMENTATION CLASS AND VERSION FOR MEDIA STORAGE**

|                             |                       |
|-----------------------------|-----------------------|
| Implementation Class UID    | 1.2.410.200001.1.0001 |
| Implementation Version Name | P6A1_0                |

## 5.2 AE SPECIFICATIONS

### 5.2.1 Offline-Media Application Entity Specification

The Offline-Media Application Entity provides standard conformance to the Media Storage Service Class. The Application Profiles and roles are listed below:

**Table 5.2-1**

**APPLICATION PROFILES, ACTIVITIES AND ROLES FOR OFFLINE-MEDIA**

| Application Profiles Supported | Real World Activity | Role     |
|--------------------------------|---------------------|----------|
| STD-US-SC-MF-DVD               | Export To DVD       | FSC, FSU |

#### 5.2.1.1 File Meta Information for the Application Entity

The File-Set Identifier included in the File Meta Header is “MEDISON P6A”.

The Source Application Entity Title included in the File Meta Header is configurable using the Setup/DICOM Menu.

#### 5.2.1.2 Real-World Activities

##### 5.2.1.2.1 Activity – Export to DVD

The Offline-Media Application Entity acts as an FSC and FSU when requested to export SOP Instances from the local database to a DVD medium.

If the contents of the current selection do not fit on a single media, a separation into multiple export jobs which can be adapted by the user will be suggested.

The user will be prompted to insert a DVD for each export job. The contents of the export job will be written together with a corresponding DICOMDIR to a DVD. Writing in multi-session mode is supported.

##### 5.2.1.2.1.1 Media Storage Application Profiles

The Offline-Media Application Entity supports the STD-US-SC-MF-DVD Application Profile.

### 5.2.1.2.1.1.1 Options

The Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the table below:

**Table 5.2-2**  
**IODS, SOP CLASSES AND TRANSFER SYNTAXES FOR OFFLINE MEDIA**

| <b>Information Object<br/>Definition</b>   | <b>SOP Class UID</b>              | <b>Transfer Syntax</b>             | <b>Transfer Syntax UID</b> |
|--|-----------------------------------|------------------------------------|----------------------------|
| Media Storage Directory<br>Storage         | 1.2.840.10008.1.3.10              | Explicit VR Little Endian          | 1.2.840.10008.1.2.1        |
| US Image Storage                           | 1.2.840.10008.5.1.4.1.1.6.1       | Explicit VR Little Endian          | 1.2.840.10008.1.2.1        |
| US Multiframe Image<br>Storage             | 1.2.840.10008.5.1.4.1.1.3.1       | JPEG Baseline Lossy<br>Compression | 1.2.840.10008.1.2.4.50     |
| Comprehensive Structured<br>Report Storage | 1.2.840.10008.5.1.4.1.1.88.3<br>3 | Explicit VR Little Endian          | 1.2.840.10008.1.2.1        |

## **6 SUPPORT OF CHARACTER SETS**

All SONOACE X8 DICOM applications support the

ISO\_IR 100 (ISO 8859-1:1987 Latin Alphabet No. 1 supplementary set)

## 7 SECURITY

SONOACE X8 does not support any specific security measures.

It is assumed that SONOACE X8 is used within a secured environment. It is assumed that a secured environment includes as minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to SONOACE X8.
- b. Firewall or router protections to ensure that SONOACE X8 has only network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriately secure network channels (e.g. such as a Virtual Private Network (VPN)).

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

## 8 ANNEXES

### 8.1 IOD CONTENTS

#### 8.1.1 Created SOP Instances

Table 8.1-1 specifies the attributes of an Ultrasound Image transmitted by the SONOACE X8 storage applications.

Table 8.1-2 specifies the attributes of a Comprehensive Structured Reports transmitted by the SONOACE X8 storage applications.

The following tables use a number of abbreviations. The abbreviations used in the “Presence of...” column are:

|        |   |
|--------|---|
| VNAP   | Value Not Always Present (attribute sends zero length if no value is present) |
| ANAP   | Attribute Not Always Present  |
| ALWAYS | Always Present  |
| EMPTY  | Attribute is sent without a value   |

The abbreviations used in the “Source” column:

|        |  |
|--------|--|
| MWL    | the attribute value source Modality Worklist                                     |
| USER   | the attribute value source is from User input                                    |
| AUTO   | the attribute value is generated automatically                                   |
| MPPS   | the attribute value is the same as the Modality Performed Procedure Step service |
| CONFIG | the attribute value source is a configurable parameter                           |

NOTE: All dates and times are encoded in the local configured calendar and time. Date, Time and Time zones are configured using the Setup Menu.

#### 8.1.1.1 US or US Multiframe Image IOD

**Table 8.1-1**  
**IOD OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| IE      | Module        | Reference   | Presence of Module |
|---------|---------------|-------------|--------------------|
| Patient | Patient       | Table 8.1-3 | ALWAYS             |
| Study   | General Study | Table 8.1-4 | ALWAYS             |

|           |                       |              |                          |
|-----------|-----------------------|--------------|--------------------------|
|           | Patient Study         | Table 8.1-5  | ALWAYS                   |
| Series    | General Series        | Table 8.1-6  | ALWAYS                   |
| Equipment | General Equipment     | Table 8.1-7  | ALWAYS                   |
| Image     | General Image         | Table 8.1-8  | ALWAYS                   |
|           | Image Pixel           | Table 8.1-9  | ALWAYS                   |
|           | Cine                  | Table 8.1-10 | Only if US<br>Multiframe |
|           | Multi-Frame           | Table 8.1-11 | Only if US<br>Multiframe |
|           | US Region Calibration | Table 8.1-12 | ANAP                     |
|           | US Image              | Table 8.1-13 | ALWAYS                   |
|           | VOI LUT               | Table 8.1-14 | ALWAYS                   |
|           | SOP Common            | Table 8.1-15 | ALWAYS                   |

### 8.1.1.2 Comprehensive Structured Report IOD

**Table 8.1-2**

#### **IOD OF CREATED COMPREHENSIVE STRUCTURED REPORT SOP INSTANCES**

| <b>IE</b> | <b>Module</b>       | <b>Reference</b> | <b>Presence of Module</b> |
|-----------|---------------------|------------------|---------------------------|
| Patient   | Patient             | Table 8.1-3      | ALWAYS                    |
| Study     | General Study       | Table 8.1-4      | ALWAYS                    |
|           | Patient Study       | Table 8.1-5      | ALWAYS                    |
| Series    | SR Document Series  | Table 8.1-16     | ALWAYS                    |
| Equipment | General Equipment   | Table 8.1-7      | ALWAYS                    |
| Document  | SR Document General | Table 8.1-17     | ALWAYS                    |
|           | SR Document Content | Table 8.1-18     | ALWAYS                    |
|           | SOP Common          | Table 8.1-19     | ALWAYS                    |

### 8.1.1.3 Common Modules

**Table 8.1-3**

#### **PATIENT MODULE OF CREATED SOP INSTANCES**

| <b>Attribute</b> | <b>Tag</b> | <b>VR</b> | <b>Value</b> | <b>Presence</b> | <b>Source</b> |
|------------------|------------|-----------|--------------|-----------------|---------------|
|------------------|------------|-----------|--------------|-----------------|---------------|

| Name                 |           |    |   | of Value |               |
|----------------------|-----------|----|---|----------|---------------|
| Patient's Name       | 0010,0010 | PN | From MWL or User Input. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain first 3 components (Last^First^Middle). Maximum 64 characters. | VNAP     | MWL/USER      |
| Patient ID           | 0010,0020 | LO | From MWL, user input or generated by device. Maximum 64 characters.   | ALWAYS   | MWL/USER/AUTO |
| Patient's Birth Date | 0010,0030 | DA | From MWL or user input  | VNAP     | MWL/USER      |
| Patient's Sex        | 0010,0040 | CS | From MWL or user input  | VNAP     | MWL/USER      |

**Table 8.1-4**

**GENERAL STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source   |
|----------------------------|-----------|----|---|-------------------|----------|
| Study Instance UID         | 0020,000D | UI | From MWL or generated by device                               | ALWAYS            | MWL/AUTO |
| Study Date                 | 0008,0020 | DA | <yyyymmdd>  | ALWAYS            | AUTO     |
| Study Time                 | 0008,0030 | TM | <hhmmss>  | ALWAYS            | AUTO     |
| Referring Physician's Name | 0008,0090 | PN | From MWL or user input  | VNAP              | MWL/USER |
| Study ID                   | 0020,0010 | SH | system generate : Study Date + Study Time<br><yyyymmddhhmmss> | ALWAYS            | AUTO     |
| Accession Number           | 0008,0050 | SH | From MWL or user input  | VNAP              | MWL/USER |
| Study Description          | 0008,1030 | LO | From MWL (Scheduled procedure step description) or user input | ANAP              | MWL/USER |

|                               |           |    |          |      |     |
|-------------------------------|-----------|----|----------|------|-----|
| Referenced Study Sequence     | 0008,1110 | SQ | From MWL | ANAP | MWL |
| > Referenced SOP Class UID    | 0008,1150 | UI | From MWL | ANAP | MWL |
| > Referenced SOP Instance UID | 0008,1155 | UI | From MWL | ANAP | MWL |
| Procedure Code Sequence       | 0008,1032 | SQ | From MWL | ANAP | MWL |

**Table 8.1-5  
PATIENT STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name   | Tag       | VR | Value                  | Presence of Value | Source   |
|------------------|-----------|----|------------------------|-------------------|----------|
| Patient's Size   | 0010,1020 | DS | From MWL or user input | ANAP              | MWL/USER |
| Patient's Weight | 0010,1030 | DS | From MWL or user input | ANAP              | MWL/USER |

**Table 8.1-6  
GENERAL SERIES MODULE OF CREATED SOP INSTANCES**

| Attribute Name      | Tag       | VR | Value               | Presence of Value | Source |
|---------------------|-----------|----|---------------------|-------------------|--------|
| Modality            | 0008,0060 | CS | US                  | ALWAYS            | AUTO   |
| Series Instance UID | 0020,000E | UI | Generated by device | ALWAYS            | AUTO   |
| Series Number       | 0020,0011 | IS | "1"                 | ALWAYS            | AUTO   |
| Series Date         | 0008,0021 | DA | <yyyymmdd>          | ALWAYS            | AUTO   |
| Series Time         | 0008,0031 | TM | <hhmmss>            | ALWAYS            | AUTO   |
| Protocol Name       | 0018,1030 | LO | "FreeForm"          | ALWAYS            | AUTO   |
| Operators' Name     | 0008,1070 | PN | From user input     | ANAP              | USER   |

|   |           |    |  |        |              |
|---|-----------|----|--|--------|--------------|
| Referenced<br>Performed<br>Procedure Step<br>Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance<br>to which this image is related | ALWAYS | MPPS         |
| > Referenced<br>SOP Class UID                         | 0008,1150 | UI | MPPS SOP Class UID<br>"1.2.840.10008.3.1.2.3.3"                    | ALWAYS | MPPS         |
| > Referenced<br>SOP Instance<br>UID                   | 0008,1155 | UI | MPPS SOP Instance UID  | ALWAYS | MPPS         |
| Request<br>Attributes<br>Sequence                     | 0040,0275 | SQ | Zero or 1 item will be present                                     | ANAP   | AUTO         |
| > Requested<br>Procedure ID                           | 0040,1001 | SH | From MWL   | ANAP   | MWL          |
| > Scheduled<br>Procedure Step<br>ID                   | 0040,0009 | SH | From MWL   | ANAP   | MWL          |
| > Scheduled<br>Procedure Step<br>Description          | 0040,0007 | LO | From MWL   | ANAP   | MWL          |
| > Scheduled<br>Protocol Code<br>Sequence              | 0040.0008 | SQ | From MWL   | ANAP   | MWL          |
| Performed<br>Procedure Step<br>ID                     | 0040,0253 | SH | Same as MPPS   | ALWAYS | MPPS         |
| Performed<br>Procedure Step<br>Start Date             | 0040,0244 | DA | Same as Study Date   | ALWAYS | AUTO         |
| Performed<br>Procedure Step<br>Start Time             | 0040,0245 | TM | Same as Study Time   | ALWAYS | AUTO         |
| Performed<br>Procedure Step<br>Description            | 0040,0254 | LO | Same as Study Description  | ANAP   | MWL/US<br>ER |

**Table 8.1-7  
GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES**

| <b>Attribute Name</b>     | <b>Tag</b> | <b>VR</b> | <b>Value</b>        | <b>Presence of Value</b> | <b>Source</b> |
|---------------------------|------------|-----------|---------------------|--------------------------|---------------|
| Manufacturer              | 0008,0070  | LO        | "MEDISON"           | ALWAYS                   | AUTO          |
| Institution Name          | 0008,0080  | LO        | From user input     | ANAP                     | CONFIG        |
| Station Name              | 0008,1010  | SH        | From user input     | ANAP                     | CONFIG        |
| Manufacturer's Model Name | 0008,1090  | LO        | "P6A"               | ALWAYS                   | AUTO          |
| Device Serial Number      | 0018,1000  | LO        | Generated by device | ALWAYS                   | AUTO          |
| Software Versions         | 0018,1020  | LO        | Generated by device | ALWAYS                   | AUTO          |

**8.1.1.4 US or US Multiframe Image Module**

**Table 8.1-8  
GENERAL IMAGE MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| <b>Attribute Name</b>         | <b>Tag</b> | <b>VR</b> | <b>Value</b>  | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------------|------------|-----------|---|--------------------------|---------------|
| Instance Number               | 0020,0013  | IS        | Generated by device, increments from "1" in each series   | ALWAYS                   | AUTO          |
| Patient Orientation           | 0020,0020  | CS        | NULL  |                          |               |
| Content Date                  | 0008,0023  | DA        | <yyyymmdd>  | ALWAYS                   | AUTO          |
| Content Time                  | 0008,0033  | TM        | <hhmmss>  | ALWAYS                   | AUTO          |
| Image Type                    | 0008,0008  | CS        | "ORIGINAL" and "PRIMARY" and "" and ""                    | ALWAYS                   | AUTO          |
| Lossy Image Compression       | 0028,2110  | CS        | US = "00" (uncompressed), US-MF = "01" (lossy compressed) | ALWAYS                   | AUTO          |
| Lossy Image Compression Ratio | 0028,2112  | DS        | Used if (0028, 2110) = "01", Calculated by device         | ANAP                     | AUTO          |
| Lossy Image Compression       | 0028,2114  | CS        | "ISO_10918_1", used if (0028,2001) = "01"                 | ANAP                     | AUTO          |

|        |  |  |  |  |  |
|--------|--|--|--|--|--|
| Method |  |  |  |  |  |
|--------|--|--|--|--|--|

**Table 8.1-9**

**IMAGE PIXEL MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name             | Tag       | VR             | Value                      | Presence of Value | Source |
|----------------------------|-----------|----------------|----------------------------|-------------------|--------|
| Samples per Pixel          | 0028,0002 | US             | "3" for RGB                | ALWAYS            | AUTO   |
| Photometric Interpretation | 0028,0004 | CS             | "RGB"                      | ALWAYS            | AUTO   |
| Rows                       | 0028,0010 | US             | US = "768", US-MF = "480"  | ALWAYS            | AUTO   |
| Columns                    | 0028,0011 | US             | US = "1024", US-MF = "640" | ALWAYS            | AUTO   |
| Bits Allocated             | 0028,0100 | US             | "8"                        | ALWAYS            | AUTO   |
| Bits Stored                | 0028,0101 | US             | "8"                        | ALWAYS            | AUTO   |
| High Bit                   | 0028,0102 | US             | "7"                        | ALWAYS            | AUTO   |
| Pixel Representation       | 0028,0103 | US             | "0"                        | ALWAYS            | AUTO   |
| Pixel Data                 | 7FE0,0010 | OW<br>or<br>OB | Generated by device        | ALWAYS            | AUTO   |
| Planar Configuration       | 0028,0006 | US             | "0"                        | ALWAYS            | AUTO   |
| Private Creator            | 7FE1,0010 | LO             | "MEDISON_US"               | ANAP              | AUTO   |
| 3D Volume                  | 7FE1,1002 | OB             | 3D Volume Data             | ANAP              | AUTO   |

**Table 8.1-10**

**CINE MODULE OF CREATED US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value             | Presence of Value | Source |
|----------------|-----------|----|-------------------|-------------------|--------|
| Frame Time     | 0018,1063 | DS | Milliseconds      | ANAP              | AUTO   |
| Cine Rate      | 0018,0040 | IS | Frames per second | ANAP              | AUTO   |

**Table 8.1-11**

**MULTI-FRAME MODULE OF CREATED US MULTIFRAME SOP INSTANCES**

| <b>Attribute Name</b>   | <b>Tag</b> | <b>VR</b> | <b>Value</b>             | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------|------------|-----------|--------------------------|--------------------------|---------------|
| Number of Frames        | 0028,0008  | IS        | Numbers of Frames        | ANAP                     | AUTO          |
| Frame Increment Pointer | 0028,0009  | AT        | "1577059" : (0018, 1063) | ANAP                     | AUTO          |

**Table 8.1-12**

**US REGION CALIBRATION MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| <b>Attribute Name</b>          | <b>Tag</b> | <b>VR</b> | <b>Value</b>   | <b>Presence of Value</b> | <b>Source</b> |
|--------------------------------|------------|-----------|--|--------------------------|---------------|
| Sequence of Ultrasound Regions | 0018,6011  | SQ        | Generated by device. A sequence is present for each region in the system display.            | ANAP                     | AUTO          |
| > Region Location Min x0       | 0018,6018  | UL        | Left position of region  | ALWAYS                   | AUTO          |
| > Region Location Min y0       | 0018,601A  | UL        | Top position of region   | ALWAYS                   | AUTO          |
| > Region Location Max x1       | 0018,601C  | UL        | Right position of region   | ALWAYS                   | AUTO          |
| > Region Location Max y1       | 0018,601E  | UL        | Bottom position of region  | ALWAYS                   | AUTO          |
| > Physical Units X Direction   | 0018,6024  | US        | 2D Image : 0003H = cm<br>M-Mode : 0004H = seconds<br>Doppler : 0004H = seconds               | ALWAYS                   | AUTO          |
| > Physical Units Y Direction   | 0018,6026  | US        | 2D Image : 0003H = cm<br>M-Mode : 0003H = cm<br>Doppler : 0005H = hertz or<br>0007H = cm/sec | ALWAYS                   | AUTO          |
| > Physical Delta X             | 0018,602C  | FD        | The physical value per pixel increment   | ALWAYS                   | AUTO          |

|                         |           |    |  |        |      |
|-------------------------|-----------|----|--|--------|------|
| > Physical Delta Y      | 0018,602E | FD | The physical value per pixel increment   | ALWAYS | AUTO |
| > Region Spatial Format | 0018,6012 | US | 2D Tissue : 0001H<br>M-Mode Tissue or flow : 0002H<br>Spectral (CW or PW Doppler) : 0003H          | ALWAYS | AUTO |
| > Region Data Type      | 0018,6014 | US | Tissue : 0001H<br>Color Flow : 0002H<br>PW Spectral Doppler : 0003H<br>CW Spectral Doppler : 0004H | ALWAYS | AUTO |
| > Region Flags          | 0018,6016 | UL | See DICOM PS 3.3 C.8.5.5.1.3   | ALWAYS | AUTO |

**Table 8.1-13**

**US IMAGE MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source |
|----------------------------|-----------|----|---|-------------------|--------|
| Samples Per Pixel          | 0028,0002 | US | "3" for RGB   | ALWAYS            | AUTO   |
| Photometric Interpretation | 0028,0004 | CS | "RGB"   | ALWAYS            | AUTO   |
| Bits Allocated             | 0028,0100 | US | "8"   | ALWAYS            | AUTO   |
| Bits Stored                | 0028,0101 | US | "8"   | ALWAYS            | AUTO   |
| High Bit                   | 0028,0102 | US | "7"   | ALWAYS            | AUTO   |
| Planar Configuration       | 0028,0006 | US | "0"   | ALWAYS            | AUTO   |
| Pixel Representation       | 0028,0103 | US | "0"   | ALWAYS            | AUTO   |
| Image Type                 | 0008,0008 | CS | "ORIGINAL" and "PRIMARY" and "" and ""                    | ALWAYS            | AUTO   |
| Lossy Image Compression    | 0028,2110 | CS | US = "00" (uncompressed), US-MF = "01" (lossy compressed) | ALWAYS            | AUTO   |

**Table 8.1-14**

**VOI LUT MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name | Tag       | VR | Value           | Presence of Value | Source |
|----------------|-----------|----|-----------------|-------------------|--------|
| Window Center  | 0028,1050 | DS | default : "128" | ALWAYS            | CONFIG |
| Window Width   | 0028,1051 | DS | default : "256" | ALWAYS            | CONFIG |

**Table 8.1-15**

**SOP COMMON MODULE OF CREATED US OR US MULTIFRAME SOP INSTANCES**

| Attribute Name         | Tag       | VR | Value   | Presence of Value | Source |
|------------------------|-----------|----|---|-------------------|--------|
| SOP Class UID          | 0008,0016 | UI | US =<br>"1.2.840.10008.5.1.4.1.1.6.1"<br>US-MF =<br>"1.2.840.10008.5.1.4.1.1.3.1" | ALWAYS            | AUTO   |
| SOP Instance UID       | 0008,0018 | UI | Generated by device   | ALWAYS            | AUTO   |
| Specific Character Set | 0008,0005 | CS | "ISO_IR 100"  | ALWAYS            | AUTO   |

**8.1.1.5 Comprehensive Structured Report Modules**

**Table 8.1-16**

**SR DOCUMENT SERIES MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| Attribute Name                               | Tag       | VR | Value   | Presence of Value | Source |
|--|-----------|----|---|-------------------|--------|
| Modality                                     | 0008,0060 | CS | SR  | ALWAYS            | AUTO   |
| Series Instance UID                          | 0020,000E | UI | Generated by device   | ALWAYS            | AUTO   |
| Series Number                                | 0020,0011 | IS | "2"   | ALWAYS            | AUTO   |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ALWAYS            | MPPS   |

|                               |           |    |   |        |      |
|-------------------------------|-----------|----|---|--------|------|
| > Referenced SOP Class UID    | 0008,1150 | UI | MPPS SOP Class UID<br>"1.2.840.10008.3.1.2.3.3" | ALWAYS | MPPS |
| > Referenced SOP Instance UID | 0008,1155 | UI | MPPS SOP Instance UID                           | ALWAYS | MPPS |

**Table 8.1-17**

**SR DOCUMENT GENERAL MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| Attribute Name                                | Tag       | VR | Value   | Presence of Value | Source   |
|---|-----------|----|---|-------------------|----------|
| Instance Number                               | 0020,0013 | IS | Generated by device, increments from "1" in each series | ALWAYS            | AUTO     |
| Completion Flag                               | 0040,A491 | CS | "PARTIAL"   | ALWAYS            | AUTO     |
| Verification Flag                             | 0040,A493 | CS | "UNVERIFIED"  | ALWAYS            | AUTO     |
| Content Date                                  | 0008,0023 | DA | <yyyymmdd>  | ALWAYS            | AUTO     |
| Content Time                                  | 0008,0033 | TM | <hhmmss>  | ALWAYS            | AUTO     |
| Referenced Request Sequence                   | 0040,A370 | SQ | 1 item will be present                                  | ALWAYS            | AUTO     |
| > Study Instance UID                          | 0020,000D | UI | From MWL or generated by device                         | ALWAYS            | MWL/AUTO |
| > Referenced Study Sequence                   | 0008,1110 | SQ | From MWL  | ANAP              | MWL      |
| >> Referenced SOP Class UID                   | 0008,1150 | UI | From MWL  | ANAP              | MWL      |
| >> Referenced SOP Instance UID                | 0008,1155 | UI | From MWL  | ANAP              | MWL      |
| > Accession Number                            | 0008,0050 | SH | From MWL or user input                                  | VNAP              | MWL/USER |
| > Placer Order Number/Imaging Service Request | 0040,2016 | LO | NULL  | VNAP              | AUTO     |

|   |           |    |          |      |      |
|---|-----------|----|----------|------|------|
| > Filler Order Number/Imaging Service Request | 0040,2017 | LO | NULL     | VNAP | AUTO |
| > Requested Procedure ID                      | 0040,1001 | SH | From MWL | VNAP | MWL  |
| > Requested Procedure Description             | 0032,1060 | LO | From MWL | VNAP | MWL  |
| > Requested Procedure Code Sequence           | 0032,1064 | SQ | From MWL | VNAP | MWL  |
| Performed Procedure Code Sequence             | 0040,A372 | SQ | NULL     | VNAP | AUTO |

**Table 8.1-18**

**SR DOCUMENT CONTENT MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| <b>Attribute Name</b>                 | <b>Tag</b> | <b>VR</b> | <b>Value</b>   | <b>Presence of Value</b> | <b>Source</b> |
|---------------------------------------|------------|-----------|--|--------------------------|---------------|
| Value Type                            | 0040,A040  | CS        | "CONTAINER"  | ALWAYS                   | AUTO          |
| Concept Name Code Sequence            | 0040,A043  | SQ        | 1 item will be present   | ALWAYS                   | AUTO          |
| > Include 'Code Sequence Macro'       |            |           | "EV(125000, DCM, "OB-GYN Ultrasound Procedure Report" for OB-GYN | ALWAYS                   | AUTO          |
| Include 'Container Macro'             |            |           |  | ALWAYS                   | AUTO          |
| Content Sequence                      | 0040,A730  | SQ        | One or more items may be included in this sequence               | ALWAYS                   | AUTO          |
| > Relationship Type                   | 0040,A010  | CS        | See ...  | ALWAYS                   | AUTO          |
| > Include Document Relationship Macro |            |           | See ...  | ALWAYS                   | AUTO          |
| > Include Document Content            |            |           | See ...  | ALWAYS                   | AUTO          |

|       |  |  |  |
|-------|--|--|--|
| Macro |  |  |  |
|-------|--|--|--|

**Table 8.1-19**

**SOP COMMON MODULE OF CREATED COMPREHENSIVE SR SOP INSTANCES**

| <b>Attribute Name</b>  | <b>Tag</b> | <b>VR</b> | <b>Value</b>                    | <b>Presence of Value</b> | <b>Source</b> |
|------------------------|------------|-----------|---------------------------------|--------------------------|---------------|
| SOP Class UID          | 0008,0016  | UI        | "1.2.840.10008.5.1.4.1.1.88.33" | ALWAYS                   | AUTO          |
| SOP Instance UID       | 0008,0018  | UI        | Generated by device             | ALWAYS                   | AUTO          |
| Specific Character Set | 0008,0005  | CS        | "ISO_IR 100"                    | ALWAYS                   | AUTO          |

**8.1.2 Used Fields in received IOD by application**

The SonoAce X8 storage application does not receive SOP Instances. The usage of attributes received via Modality Worklist is described in section 4.2.2.3.1.3.

**8.1.3 Attribute mapping**

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in the Table below. The format and conversions used in Table are the same as the corresponding table in IHE Technical Framework, Rev. 7.0 May 15, 2006, vol. II, Appendix A.

**Table 8.1-20**

**ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST, IMAGE AND MPPS**

| <b>Modality Worklist</b>   | <b>Image IOD</b>           | <b>MPPS IOD</b>                    |
|----------------------------|----------------------------|------------------------------------|
| Patient's Name             | Patient's Name             | Patient's Name                     |
| Patient ID                 | Patient ID                 | Patient ID                         |
| Patient's Birth Date       | Patient's Birth Date       | Patient's Birth Date               |
| Patient's Sex              | Patient's Sex              | Patient's Sex                      |
| Patient's Size             | Patient's Size             | _____                              |
| Patient's Weight           | Patient's Weight           | _____                              |
| Referring Physician's Name | Referring Physician's Name | _____                              |
|                            | _____                      | Scheduled Step Attributes Sequence |

|                                      |  |  |
|--------------------------------------|--|--|
| Study Instance UID                   | Study Instance UID                           | > Study Instance UID                   |
| Referenced Study Sequence            | Referenced Study Sequence                    | > Referenced Study Sequence            |
| Accession Number                     | Accession Number                             | > Accession Number                     |
|                                      | Request Attributes Sequence                  | _____                                  |
| Requested Procedure ID               | > Requested Procedure ID                     | > Requested Procedure ID               |
| Requested Procedure Description      | _____  | > Requested Procedure Description      |
| Scheduled Procedure Step ID          | > Scheduled Procedure Step ID                | > Scheduled Procedure Step ID          |
| Scheduled Procedure Step Description | > Scheduled Procedure Step Description       | > Scheduled Procedure Step Description |
| Scheduled Protocol Code Sequence     | > Scheduled Protocol Code Sequence           | > Scheduled Protocol Code Sequence     |
|                                      |  |  |
| _____                                | Study ID                                     | Study ID                               |
| _____                                | Performed Procedure Step ID                  | Performed Procedure Step ID            |
| _____                                | Performed Procedure Step Start Date          | Performed Procedure Step Start Date    |
| _____                                | Performed Procedure Step Start Time          | Performed Procedure Step Start Time    |
| _____                                | Performed Procedure Step Description         | Performed Procedure Step Description   |
| _____                                | _____  | Performed Series Sequence              |
| Requested Procedure Code Sequence    | Procedure Code Sequence                      | Procedure Code Sequence                |
| _____                                | Referenced Performed Procedure Step Sequence | _____                                  |
| _____                                | > Referenced SOP Class UID                   | SOP Class UID                          |
| _____                                | > Referenced SOP Instance UID                | SOP Instance UID                       |
| _____                                | Protocol Name                                | Protocol Name                          |

#### 8.1.4 Coerced/Modified Fields

The Modality Worklist AE will truncate attribute values received in the response to a Modality Worklist Query if the value length is longer than the maximum length permitted by the attribute's VR.

## 8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

The Private Attributes added to create SOP Instances are listed in the Table below. SONOACE X8 reserves blocks of private attributes in groups 7FE1. Further details on usage of these private attributes are contained in Section 8.1

**Table 8.2-1**  
**DATA DICTIONARY OF PRIVATE ATTRIBUTES**

| Tag          | Attribute Name  | VR | VM |
|--------------|-----------------|----|----|
| (7FE1, 0010) | Private Creator | LO | 1  |
| (7FE1, 1002) | 3D Volume       | OB | 1  |

## 8.3 CODED TERMINOLOGY AND TEMPLATES

The Workflow AE is capable of supporting arbitrary coding schemes for Procedure and Protocol Codes. The contents of Requested Procedure Code Sequence (0032, 1064) and Scheduled Protocol Code Sequence (0040, 0008) supplied in Worklist Items will be mapped to Image IOD and MPPS attributes as described in Section 8.1.3

## 8.4 GRAYSCALE IMAGE CONSISTENCY

The high resolution display monitor attached to SONOACE X8 can be calibrated according to the Grayscale Standard Display Function (GSDF).

## 8.5 STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES

No Specialized or Private SOP Classes are supported.

### 8.5.1 US OR US MULTIFRAME IMAGE STORAGE SOP CLASS

The US or US Multiframe Image Storage SOP Classes are extended to create a Standard Extended SOP Class by addition of standard and private attributes to the created SOP Instances as documented in section 8.1

3D Volume Data is transferred to the configured Storage Server, if "Send 3D Volume" option is enabled in the Setup Dialog.

## 8.6 PRIVATE TRANSFER SYNTAXES

No Private Transfer Syntaxes are supported.

## 9 STRUCTURED REPORT TEMPLATES

### 9.1 TEMPLATES used in SONOACE X8

This Section uses the following forms for describing Structured Report Templates used in SONOACE X8.

|   | Rel with Parent | VT | Concept Name | Presence of Value | Comments |
|---|-----------------|----|--------------|-------------------|----------|
| 1 |                 |    |              |                   |          |
| 2 |                 |    |              |                   |          |

|     | NL | REL | VT | Concept Name | Unit / CODE Value | Ref TID | Ref CID | Comments |
|-----|----|-----|----|--------------|-------------------|---------|---------|----------|
| A-1 |    |     |    |              |                   |         |         |          |
| A-2 |    |     |    |              |                   |         |         |          |

Rel with Parent      Relationship

VT                      Value Type

Concept Name      Any constraints on Concept Name are specified in this field as defined or enumerated coded entries, or as baseline or defined context groups.

Presence of Value   Ref. Section 8.1.1

Comments            Description about Reference section or used values.

NL                      The nesting level of Content Items is denoted by ">" symbols

REL                    Relationship

Unit/Code, Value    Applied unit, enumerated coded entries, or the reference of Context Group.

Ref TID                Referenced Template ID Number

Ref CID                Referenced Context ID Number. The left side of "/" shows a CID value applied in "Concept Name" column and the right side shows a CID value applied in "Unit/Code, Value" column. (e.g. 228/12012)

## 9.1.1 OB-GYN STRUCTURED REPORT TEMPLATE


### 9.1.1.1 OB-GYN Ultrasound Procedure Report (TID 5000)

**Table 9.1-1  
OB-GYN ULTRASOUND PROCEDURE REPORT TEMPLATE**

|    | Rel with Parent | VT        | Concept Name   | Presence of Value | Comments              |
|----|-----------------|-----------|--|-------------------|-----------------------|
| 1  |                 | CONTAINER | EV (125000, DCM, "OB-GYN Ultrasound Procedure Report") | ALWAYS            |                       |
| 2  | HAS CONCEPT MOD | INCLUDE   | DTID (1204) Language of Content Item and Descendants   |                   |                       |
| 3  | HAS OBS CONTEXT | INCLUDE   | DTID (1001) Observation Context                        | ANAP              | Ref. Section 9.1.1.2  |
| 4  | CONTAINS        | INCLUDE   | DTID (5001) Patient Characteristics                    | ANAP              | Ref. Section 9.1.1.3  |
| 5  | CONTAINS        | CONTAINER | DT (111028, DCM, "Image Library")                      |                   |                       |
| 6  | CONTAINS        | IMAGE     | No Purpose of reference                                |                   |                       |
| 7  | CONTAINS        | INCLUDE   | DTID (5002) OB-GYN Procedure Summary Section           | ANAP              | Ref. Section 9.1.1.4  |
| 8  | CONTAINS        | INCLUDE   | DTID (5004) Fetal Biometry Ratio Section               | ANAP              | Ref. Section 9.1.1.5  |
| 9  | CONTAINS        | INCLUDE   | DTID (5005) Fetal Biometry Section                     | ANAP              | Ref. Section 9.1.1.6  |
| 10 | CONTAINS        | INCLUDE   | DTID (5006) Long Bones Section                         | ANAP              | Ref. Section 9.1.1.7  |
| 11 | CONTAINS        | INCLUDE   | DTID (5007) Fetal Cranium Section                      | ANAP              | Ref. Section 9.1.1.8  |
| 12 | CONTAINS        | INCLUDE   | DTID (5009) Fetal Biophysical Profile Section          | ANAP              | Ref. Section 9.1.1.9  |
| 13 | CONTAINS        | INCLUDE   | DTID (5011) Early Gestation Section                    | ANAP              | Ref. Section 9.1.1.11 |
| 14 | CONTAINS        | INCLUDE   | DTID (5010) Amniotic Sac                               | ANAP              | Ref. Section 9.1.1.10 |

|    |                 |           |  |      |                       |
|----|-----------------|-----------|--|------|-----------------------|
|    |                 |           | Section  |      |                       |
| 15 | CONTAINS        | INCLUDE   | DTID (5015) Pelvis and Uterus Section                | ANAP | Ref. Section 9.1.1.12 |
| 16 | CONTAINS        | INCLUDE   | DTID (5012) Ovaries Section                          | ANAP | Ref. Section 9.1.1.13 |
| 17 | CONTAINS        | INCLUDE   | DTID (5013) Follicles Section                        | ANAP | Ref. Section 9.1.1.14 |
| 18 | CONTAINS        | INCLUDE   | DTID (5013) Follicles Section                        | ANAP | Ref. Section 9.1.1.15 |
| 19 | CONTAINS        | CONTAINER | EV (121070, DCM, "Findings")                         |      |                       |
| 20 | HAS CONCEPT MOD | CODE      | EV (G-C0E3, SRT, "Finding Site")                     |      |                       |
| 21 | CONTAINS        | INCLUDE   | DTID (5025) OB-GYN Fetal Vascular Measurement Group  |      |                       |
| 22 | CONTAINS        | CONTAINER | EV (121070, DCM, "Findings")                         |      |                       |
| 23 | HAS CONCEPT MOD | CODE      | EV (G-C0E3, SRT, "Finding Site")                     |      |                       |
| 24 | CONTAINS        | INCLUDE   | DTID (5026) OB-GYN Pelvic Vascular Measurement Group |      |                       |

### 9.1.1.2 Observation Context (TID 1001)




**Figure 9.1-1**  
**TEMPLATE HIERARCHY OF OBSERVATION CONTEXT IN OB-GYN SR**

**Table 9.1-2  
OBSERVATION CONTEXT IN OB-GYN SR**

|     | NL | REL             | VT    | Concept Name                       | Unit / CODE Value                                  | Ref TID | Ref CID | Comments  |
|-----|----|-----------------|-------|------------------------------------|--|---------|---------|---|
| A-1 | >  | HAS OBS CONTEXT | CODE  | DCM 121005<br>Observer Type        | DCM 121006<br>"Person"                             | 1002    | /270    | This will have a value "Person".  |
| A-2 | >  | HAS OBS CONTEXT | PNAME | DCM 121008<br>Person Observer Name |  | 1003    |         | Value is taken from "Ref. Physician" of Study Info dialog.                  |
| A-3 | >  | HAS OBS CONTEXT | CODE  | DCM 121024<br>Subject Class        | DCM 121025<br>"Patient"                            | 1006    | /271    | This will have a value "Patient".   |
| A-4 | >  | HAS OBS CONTEXT | PNAME | DCM 121029<br>Subject Name         |  | 1007    |         | Value is taken from "Last Name" and "First Name" of Patient Manager dialog. |
| A-5 | >  | HAS OBS CONTEXT | DATE  | DCM 121031<br>Subject Birth Date   | yyyymmdd   |         |         | Value is taken from "Birth" of Patient Manager dialog.                      |
| A-6 | >  | HAS OBS CONTEXT | CODE  | DCM 121032<br>Subject Sex          | DCM M Male<br>DCM F Female<br>DCM U<br>Unknown sex |         |         | /7455   |
| A-7 | >  | HAS OBS CONTEXT | NUM   | DCM 121033<br>Subject Age          | UCUM mo<br>month                                   |         | 7456    | Not used  |

**9.1.1.3 Patient Characteristics (TID 5001)**




**Figure 9.1-2  
TEMPLATE HIERARCHY OF PATIENT CHARACTERISTICS IN OB-GYN SR**

**Table 9.1-3  
PATIENT CHARACTERISTICS IN OB-GYN SR**

|       | NL | REL      | VT        | Concept Name                             | Unit / CODE<br>Value                           | Ref<br>TID | Ref<br>CID | Comments  |
|-------|----|----------|-----------|--|--|------------|------------|---|
| A-8   | >  | CONTAINS | CONTAINER | DCM 121118<br>Patient<br>Characteristics |  | 5001       |            |   |
| A-8-1 | >> | CONTAINS | TEXT      | DCM 121106<br>Comment                    |  |            |            | Value is taken from<br>"Description" of<br>Study Info dialog. |
| A-8-2 | >> | CONTAINS | NUM       | LN 8302-2 Patient<br>Height              | UCUM cm<br>centimeter<br>UCUM mm<br>millimeter |            |            | Value is taken from<br>Patient Manager<br>dialog.             |
| A-8-3 | >> | CONTAINS | NUM       | LN 29463-7<br>Patient Weight             | UCUM kg<br>kilograms                           |            |            | Value is taken from<br>Patient Manager<br>dialog.             |
| A-8-4 | >> | CONTAINS | NUM       | LN 11996-6<br>Gravida                    | UCUM 1 no<br>units                             |            |            | Value is taken from<br>Category OB of<br>Study Info dialog.   |
| A-8-5 | >> | CONTAINS | NUM       | LN 11977-6 Para                          | UCUM 1 no<br>units                             |            |            | Value is taken from<br>Category OB of<br>Study Info dialog.   |
| A-8-6 | >> | CONTAINS | NUM       | LN 11612-9<br>Aborta                     | UCUM 1 no<br>units                             |            |            | Value is taken from<br>Category OB of<br>Study Info dialog.   |
| A-8-7 | >> | CONTAINS | NUM       | LN 33065-4<br>Ectopic<br>Pregnancies     | UCUM 1 no<br>units                             |            |            | Value is taken from<br>Category OB of<br>Study Info dialog.   |

9.1.1.4 OB-GYN Procedure Summary Section (TID 5002)



**Figure 9.1-3**  
**TEMPLATE HIERARCHY OF OB-GYN PROCEDURE SUMMARY SECTION**

**Table 9.1-4**  
**OB-GYN PROCEDURE SUMMARY SECTION**

|        | NL | REL      | VT        | Concept Name               | Unit / CODE<br>Value | Ref<br>TID | Ref<br>CID | Comments   |
|--------|----|----------|-----------|----------------------------|----------------------|------------|------------|--|
| A-10   | >  | CONTAINS | CONTAINER | DCM 121111<br>Summary      |                      | 5002       |            |  |
| A-10-1 | >> | CONTAINS | DATE      | LN 11778-8 EDD             | yyyymmdd             |            | 12003      | Value is taken from<br>Category OB of<br>Study Info dialog.<br>This Name for X8 is<br>"Estab.DueDate". |
|        |    |          |           | LN 11779-6 EDD<br>from LMP | yyyymmdd             |            |            | Value automatically<br>calculated by the   |


X8 based on the  
value entered for  
LMP.

|        |    |          |      |  |                    |      |       |  |
|--------|----|----------|------|--|--------------------|------|-------|--|
|        |    |          |      |  |                    |      |       | This Name for X8 is "EDD(LMP)"   |
|        |    |          |      | LN 11781-2 EDD<br>from average<br>ultrasound age   | yyymmdd            |      |       | Value automatically calculated by the X8 based various measurements and on the LMP. If there is more than one fetus, the value used is EDD of Fetus A.<br>This Name for X8 is "EDD(Average US GA)" |
|        |    |          |      | LN 11955-2 LMP                                     | yyymmdd            |      |       | Value is taken from Category OB of Study Info dialog.  |
|        |    |          |      | LN 11976-8<br>Ovulation date                       | yyymmdd            |      |       | Value is taken from Category OB of Study Info dialog.<br>This Name for X8 is "Exp.Ovul."   |
| A-10-2 | >> | CONTAINS | NUM  | LN 11878-6<br>Number of<br>Fetuses                 | UCUM 1 no<br>units | 300  | 12018 | Value is taken from Category OB of Study Info dialog.<br>This Name for X8 is "Gestations".   |
|        |    |          |      | LN 11886-9<br>Gestational Age<br>by ovulation date |                    |      |       | Not used   |
| A-10-3 | >> | CONTAINS | TEXT | DCM 121106<br>Comment                              |                    | 5002 |       | This field is taken from "Comment" entered in the  |

|            |      |                    |           |   |   |      |                          |   |
|------------|------|--------------------|-----------|---|---|------|--------------------------|---|
|            |      |                    |           |   |   |      |                          | Report.   |
| A-10-4     | >>   | CONTAINS           | CONTAINER | DCM 125008<br>Fetus Summary                             |   | 5003 |                          | This template is included 1 per fetus.<br>X8 used this template to insert measurements from DCID 12019. |
| A-10-4-1   | >>>  | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID                                  |   | 1008 |                          | Value of "1, "2, "3 or "4 is used as identifier of the Fetus.   |
| A-10-4-2   | >>>  | CONTAINS           | NUM       | LN 11878-6<br>Number of<br>Fetuses                      |   |      |                          | Not used  |
| A-10-4-3   | >>>  | CONTAINS           | NUM       | LN 18185-9<br>Gestational Age                           | UCUM d days   | 300  | 12019                    | This is a system-calculated value.<br>This name for X8 is "Average US GA"                               |
|            |      |                    |           | LN 11885-1<br>Gestational Age<br>by LMP                 | UCUM d days   |      |                          | Value automatically calculated by the X8 based on the value entered for LMP.                            |
| A-10-4-4   | >>>  | CONTAINS           | NUM       | LN 11727-5<br>Estimated Weight                          | UCUM kg   |      |                          | This is a system-calculated value.<br>This name for X8 is "EFW"   |
| A-10-4-4-1 | >>>> | INFERRED<br>FROM   | CODE      | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.2 OB Fetal<br>Body Weight<br>Equations and<br>Tables<br>(Context ID | 315  | 228<br>/12012(<br>12014) | CID 12014 will be used.   |

|            |      |               |      |   |   |     |                   |  |
|------------|------|---------------|------|---|---|-----|-------------------|--|
|            |      |               |      |   | 12014)  |     |                   |  |
| A-10-4-5   | >>>  | CONTAINS      | NUM  | LN 11767-1 EFW percentile rank                    | UCUM percentile "percentile"  | 300 | 12019             | This is a system-calculated value. This name for x8 is "Percentile (EFW)". |
| A-10-4-5-1 | >>>> | INFERRED FROM | CODE | DCM 121420 Equation<br>DCM 121424 Table of Values | Ref. Section 9.2.4<br>Estimated Fetal Weight Percentile Equations and Tables (Context ID 12016) | 315 | 228 /12012(12016) | CID 12016 will be used.  |
| A-10-4-6   | >>>  | CONTAINS      | NUM  | LN 11948-7 Fetal Heart Rate                       | UCUM bpm "bpm"  | 300 | 12019             | This is a measured value. This name for X8 is "FHR"                        |

### 9.1.1.5 Fetal Biometry Ratio Section (TID 5004)




**Figure 9.1-4**

### TEMPLATE HIERARCHY OF FETAL BIOMETRY RATIO SECTION IN OB-GYN SR

**Table 9.1-5  
FETAL BIOMETRY RATIO SECTION IN OB-GYN SR**

|        | NL | REL                | VT        | Concept Name                           | Unit / CODE<br>Value | Ref<br>TID | Ref<br>CID | Comments  |
|--------|----|--------------------|-----------|--|----------------------|------------|------------|---|
| A-11   | >  | CONTAINS           | CONTAINER | DCM 125001<br>Fetal Biometry<br>Ratios |                      | 5004       |            | Measurements<br>from CID 12004<br>are included.                         |
| A-11-1 | >> | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID                 |                      | 1008       |            | Value of "1, "2, "3,<br>or "4 is used as<br>identifier of the<br>Fetus. |
| A-11-2 | >> | CONTAINS           | NUM       | LN 11947-9<br>HC/AC                    | UCUM 1 no<br>units   | 5004       | 12004      | HC/AC   |
|        |    |                    |           | LN 11871-1<br>FL/AC                    | UCUM % "%"           |            |            | FL/AC   |
|        |    |                    |           | LN 11872-9<br>FL/BPD                   | UCUM % "%"           |            |            | FL/BPD  |
|        |    |                    |           | LN 11823-2<br>Cephalic Index           | UCUM % "%"           |            |            | CI(BPD/OFD)   |
|        |    |                    |           | LN 11873-7<br>FL/HC                    | UCUM % "%"           |            |            | FL/HC   |

**9.1.1.6 Fetal Biometry Section (TID 5005)**



**Figure 9.1-5**  
**TEMPLATE HIERARCHY OF FETAL BIOMETRY SECTION IN OB-GYN SR**


**Table 9.1-6**  
**FETAL BIOMETRY SECTION IN OB-GYN SR**

|          | NL  | REL                | VT        | Concept Name  | Unit / CODE<br>Value  | Ref<br>TID | Ref<br>CID | Comments  |
|----------|-----|--------------------|-----------|---|-----------------------|------------|------------|---|
| A-12     | >   | CONTAINS           | CONTAINER | DCM 125002<br>Fetal Biometry                              |                       | 5005       |            |   |
| A-12-1   | >>  | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID                                    |                       | 1008       |            | Will be present if<br>more than one<br>fetus.   |
| A-12-2   | >>  | CONTAINS           | CONTAINER | DCM 125005<br>Biometry Group                              |                       | 5008       |            | Measurements<br>from DCID 12005<br>are used to invoke<br>this template one<br>or more number of<br>times. |
| A-12-2-1 | >>> | CONTAINS           | NUM       | LN 11979-2<br>Abdominal<br>Circumference                  | UCUM cm<br>centimeter | 300        | 12005      | AC  |
|          |     |                    |           | LN 11818-2<br>Anterior-Posterior<br>Abdominal<br>Diameter | UCUM cm<br>centimeter |            |            | APD   |

|  |  |  |  |  |  |            |
|--|--|--|--|--|--|------------|
|  |  |  | LN 11819-0<br>Anterior-Posterior<br>Trunk Diameter | UCUM cm<br>centimeter                        |  | APTD       |
|  |  |  | LN 11820-8<br>Biparietal<br>Diameter               | UCUM cm<br>centimeter                        |  | BPD        |
|  |  |  | LN 11963-6<br>Femur Length                         | UCUM cm<br>centimeter                        |  | FL         |
|  |  |  | LN 11965-1 Foot<br>length                          | UCUM cm<br>centimeter                        |  | Foot       |
|  |  |  | LN 11984-2 Head<br>Circumference                   | UCUM cm<br>centimeter                        |  | HC         |
|  |  |  | LN 11851-3<br>Occipital-Frontal<br>Diameter        | UCUM cm<br>centimeter                        |  | OFD        |
|  |  |  | LN 11988-3<br>Thoracic<br>Circumference            | UCUM cm<br>centimeter                        |  | ThC        |
|  |  |  | LN 33068-8<br>Thoracic Area                        | UCUM cm <sup>2</sup><br>Square<br>centimeter |  | ThA        |
|  |  |  | LN 11862-0<br>Transverse<br>Abdominal<br>Diameter  | UCUM cm<br>centimeter                        |  | TAD        |
|  |  |  | LN 11864-6<br>Transverse<br>Thoracic<br>Diameter   | UCUM cm<br>centimeter                        |  | TTD        |
|  |  |  | LN 11834-9 Left<br>Kidney length                   | UCUM cm<br>centimeter                        |  | Lt. Kidney |
|  |  |  | LN 11836-4 Right<br>Kidney length                  | UCUM cm<br>centimeter                        |  | Rt. Kidney |

|            |      |                       |      |   |   |               |  |
|------------|------|-----------------------|------|---|---|---------------|--|
| A-12-2-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation                                | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean  | /3627         | If user selects Avg.<br>from the Report,<br>this value will be a<br>"Mean". Else this<br>value will be a<br>"Best value" |
| A-12-2-2   | >>>  | CONTAINS              | NUM  | LN 18185-9<br>Gestational Age                           | UCUM d days   |               | This value<br>automatically<br>calculates the GA<br>based on GA<br>equations and GA<br>tables.                           |
| A-12-2-2-1 | >>>> | INFERRED<br>FROM      | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.1<br>Gestational<br>Age Equations<br>and Tables<br>(Context<br>Group 12013) | 228<br>/12013 | CID 12013 will be<br>used.   |
| A-12-2-3   | >>>  | CONTAINS              | NUM  | DCM 125012<br>Growth Percentile<br>Rank                 | UCUM<br>percentile<br>"percentile"  | 5008<br>12017 | This value<br>automatically<br>calculates the FG<br>percentile based<br>on FG equations<br>and FG tables.                |
| A-12-2-3-1 | >>>> | INFERRED<br>FROM      | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.3 Fetal<br>Growth<br>Equations and<br>Tables<br>(Context ID<br>12015)       | 228<br>/12015 | CID 12015 will be<br>used.   |

### 9.1.1.7 Long Bones Section (TID 5006)



**Figure 9.1-6**  
**TEMPLATE HIERARCHY OF LONG BONES SECTION IN OB-GYN SR**


**Table 9.1-7**  
**LONG BONES SECTION IN OB-GYN SR**

|          | NL  | REL                | VT        | Concept Name                   | Unit / CODE Value     | Ref TID | Ref CID | Comments  |
|----------|-----|--------------------|-----------|--------------------------------|-----------------------|---------|---------|---|
| A-13     | >   | CONTAINS           | CONTAINER | DCM 125003<br>Fetal Long Bones |                       | 5006    |         |   |
| A-13-1   | >>  | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID         |                       | 1008    |         | Will be present if<br>more than one<br>fetus.   |
| A-13-2   | >>  | CONTAINS           | CONTAINER | DCM 125005<br>Biometry Group   |                       | 5008    |         | Measurements<br>from DCID 12006<br>are used to invoke<br>this template one<br>or more number of<br>times. |
| A-13-2-1 | >>> | CONTAINS           | NUM       | LN 11966-9<br>Humerus length   | UCUM cm<br>centimeter | 300     | 12006   | HUM   |
|          |     |                    |           | LN 11967-7<br>Radius length    |                       |         |         | RAD   |

|            |      |                       |      |   |   |      |        |   |
|------------|------|-----------------------|------|---|---|------|--------|---|
|            |      |                       |      | LN 11969-3 Ulna length                                  |   |      |        | ULNA  |
|            |      |                       |      | LN 11968-5 Tibia length                                 |   |      |        | TIB   |
|            |      |                       |      | LN 11964-4 Fibula length                                |   |      |        | FIB   |
|            |      |                       |      | LN 11962-8 Clavicle length                              |   |      |        | CLAV  |
| A-13-2-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation                                | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean  | 300  | /3627  | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value" |
| A-13-2-2   | >>>  | CONTAINS              | NUM  | LN 18185-9<br>Gestational Age                           | UCUM d day  | 5008 |        | This value is automatically calculates the GA based on GA equations and GA tables.                        |
| A-13-2-2-1 | >>>> | INFERRED<br>FROM      | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.1<br>Gestational<br>Age Equations<br>and Tables<br>(Context<br>Group 12013) | 228  | /12013 | CID 12013 will be used.   |
| A-13-2-3   | >>>  | CONTAINS              | NUM  | DCM 125012<br>Growth Percentile<br>Rank                 | UCUM<br>percentile<br>"percentile"  |      | 12017  | This value automatically calculates the FG percentile based on FG equations and FG tables.                |

|            |      |               |      |   |   |               |                         |
|------------|------|---------------|------|---|---|---------------|-------------------------|
| A-13-2-3-1 | >>>> | INFERRED FROM | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.3 Fetal<br>Growth<br>Equations and<br>Tables<br>(Context ID<br>12015) | 228<br>/12015 | CID 12015 will be used. |
|------------|------|---------------|------|---|---|---------------|-------------------------|

### 9.1.1.8 Fetal Cranium Section (TID 5007)



**Figure 9.1-7**  
**TEMPLATE HIERARCHY OF FETAL CRANIUM SECTION IN OB-GYN SR**


**Table 9.1-8**  
**FETAL CRANIUM SECTION IN OB-GYN SR**

|        | NL | REL                | VT        | Concept Name                | Unit / CODE<br>Value | Ref<br>TID | Ref<br>CID | Comments                                |
|--------|----|--------------------|-----------|-----------------------------|----------------------|------------|------------|---|
| A-14   | >  | CONTAINS           | CONTAINER | DCM 125004<br>Fetal Cranium |                      | 5007       |            |   |
| A-14-1 | >> | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID      |                      | 1008       |            | Will be present if more than one fetus. |

|            |      |                       |           |   |   |       |   |
|------------|------|-----------------------|-----------|---|---|-------|---|
| A-14-2     | >>   | CONTAINS              | CONTAINER | DCM 125005<br>Biometry Group  |   | 5008  | Measurements from DCID 12007 are used to invoke this template one or more of times.                       |
| A-14-2-1   | >>>  | CONTAINS              | NUM       | LN 12171-5<br>Lateral Ventrical width<br>LN 11860-4<br>Cisterna Magna length<br>LN 12146-7<br>Nuchal Fold thickness<br>LN 33070-4 Inner<br>Orbital Diameter<br>LN 11629-3 Outer<br>Orbital Diameter<br>LN 11863-8 Trans<br>Cerebellar Diameter<br>LN 33069-6<br>Nuchal Translucency | UCUM cm<br>centimeter<br>UCUM cm2<br>Square<br>centimeter | 300   | 12007<br>Lat Vent<br>CM<br>NF<br>IOD<br>OOD<br>CEREB<br>Nuchal Thickness                                  |
| A-14-2-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE      | DCM 121401<br>Derivation  | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean          | /3627 | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value" |
| A-14-2-2   | >>>  | CONTAINS              | NUM       | LN 18185-9<br>Gestational Age   | UCUM d day  | 5008  | This value automatically calculates the GA based on GA equations and GA tables.                           |

|            |      |               |      |   |   |               |  |
|------------|------|---------------|------|---|---|---------------|--|
| A-14-2-2-1 | >>>> | INFERRED FROM | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.1<br>Gestational<br>Age Equations<br>and Tables<br>(Context<br>Group 12013) | 228<br>/12013 | CID 12013 will be used.  |
| A-14-2-3   | >>>  | CONTAINS      | NUM  | DCM 125012<br>Growth Percentile<br>Rank                 | UCUM<br>percentile<br>"percentile"  | 12017         | This value automatically calculates the FG percentile based on FG equations and FG tables. |
| A-14-2-3-1 | >>>> | INFERRED FROM | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.3 Fetal<br>Growth<br>Equations and<br>Tables<br>(Context ID<br>12015)       | 228<br>/12015 | CID 12015 will be used.  |

### 9.1.1.9 Fetal Biophysical Profile Section (TID 5009)




**Figure 9.1-8**

### TEMPLATE HIERARCHY OF FETAL BIOPHYSICAL PROFILE SECTION IN OB-GYN SR

**Table 9.1-9  
FETAL BIOPHYSICAL PROFILE SECTION IN OB-GYN SR**

|        | NL | REL                | VT        | Concept Name                                   | Unit / CODE<br>Value         | Ref<br>TID | Ref<br>CID | Comments   |
|--------|----|--------------------|-----------|--|------------------------------|------------|------------|--|
| A-16   | >  | CONTAINS           | CONTAINER | DCM 125006<br>Biophysical<br>Profile           |                              | 5009       |            |  |
| A-16-1 | >> | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID                         |                              | 1008       |            | Will be present if<br>more than one<br>fetus.            |
| A-16-2 | >> | CONTAINS           | NUM       | LN 11631-9<br>Gross Body<br>Movement           | UCUM<br>{0:2} "range<br>0:2" | 5009       |            | X8 uses the value<br>as entered in the<br>Report.        |
|        |    |                    |           | LN 11632-7 Fetal<br>Breathing                  |                              |            |            | X8 uses the value<br>as entered in the<br>Report.        |
|        |    |                    |           | LN 11635-0 Fetal<br>Tone                       |                              |            |            | X8 uses the value<br>as entered in the<br>Report.        |
|        |    |                    |           | LN 11635-5 Fetal<br>Heart Reactivity           |                              |            |            | X8 uses the value<br>as entered in the<br>Report.        |
|        |    |                    |           | LN 11630-1<br>Amniotic Fluid<br>Volume         |                              |            |            | X8 uses the value<br>as entered in the<br>Report.        |
|        |    |                    |           | LN 11634-3<br>Biophysical<br>Profile Sum Score | UCUM 1 no<br>units           |            |            | X8 automatically<br>calculates the sum<br>of the scores. |

**9.1.1.10 Early Gestation Section (TID 5011)**



**Figure 9.1-9**

**TEMPLATE HIERARCHY OF EARLY GESTATION SECTION IN OB-GYN SR**


**Table 9.1-10**

**EARLY GESTATION SECTION IN OB-GYN SR**

|          | NL  | REL                | VT        | Concept Name                              | Unit / CODE<br>Value  | Ref<br>TID | Ref<br>CID | Comments  |
|----------|-----|--------------------|-----------|---|-----------------------|------------|------------|---|
| A-15     | >   | CONTAINS           | CONTAINER | DCM, 125009<br>Early Gestation            |                       | 5011       |            |   |
| A-15-1   | >>  | HAS OBS<br>CONTEXT | TEXT      | LN 11951-1 Fetus<br>ID                    |                       | 1008       |            | Will be present if<br>more than one<br>fetus.   |
| A-15-2   | >>  | CONTAINS           | CONTAINER | DCM 125005<br>Biometry Group              |                       | 5008       |            | Measurements<br>from DCID 12009<br>are used to invoke<br>this template one<br>or more number of<br>times. |
| A-15-2-1 | >>> | CONTAINS           | NUM       | LN 11957-8<br>Crown Rump<br>Length        | UCUM cm<br>centimeter | 300        | 12009      | CRL   |
|          |     |                    |           | LN 11850-5<br>Gestational Sac<br>Diameter | Square<br>centimeter  |            |            | GS  |
|          |     |                    |           | LN 33071-2 Spine<br>Length                |                       |            |            | SL  |

|            |      |                       |      |   |   |               |   |
|------------|------|-----------------------|------|---|---|---------------|---|
|            |      |                       |      | LN 11816-6 Yolk<br>Sac length                           |   |               | YS  |
| A-15-2-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation                                | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean  | /3627         | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value" |
| A-15-2-2   | >>>  | CONTAINS              | NUM  | LN 18185-9<br>Gestational Age                           | UCUM d day  |               | This value is automatically calculates the GA based on GA equations and GA tables.                        |
| A-15-2-2-1 | >>>> | INFERRED<br>FROM      | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.1<br>Gestational<br>Age Equations<br>and Tables<br>(Context Group<br>12013) | 228<br>/12013 | CID 12013 will be used.   |
| A-15-2-3   | >>>  | CONTAINS              | NUM  | DCM 125012<br>Growth Percentile<br>Rank                 | UCUM<br>percentile<br>"percentile"  | 12017         | This value automatically calculates the FG percentile based on FG equations and FG tables.                |
| A-15-2-3-1 | >>>> | INFERRED<br>FROM      | CODE | DCM 121420<br>Equation<br>DCM 121424<br>Table of Values | Ref. Section<br>9.2.3 Fetal<br>Growth<br>Equations and<br>Tables<br>(Context ID<br>12015)       | 228<br>/12015 | CID 12015 will be used.   |

9.1.1.11 Amniotic Sac Section (TID 5010)




**Figure 9.1-10**  
**TEMPLATE HIERARCHY OF AMNIOTIC SAC SECTION IN OB-GYN SR**

**Table 9.1-11**  
**AMNIOTIC SAC SECTION IN OB-GYN SR**

|        | NL | REL                   | VT        | Concept Name                              | Unit / CODE<br>Value              | Ref<br>TID | Ref<br>CID | Comments                             |
|--------|----|-----------------------|-----------|---|-----------------------------------|------------|------------|--------------------------------------|
| A-17   | >  | CONTAINS              | CONTAINER | DCM 121070<br>Findings                    |                                   |            |            |                                      |
| A-17-1 | >> | HAS<br>CONCEPT<br>MOD | CODE      | SRT G-C0E3<br>Finding Site                | SRT T-<br>F1300 "Amniotic<br>Sac" | 5010       |            | DT (T-F1300, SRT,<br>"Amniotic Sac") |
| A-17-2 | >> | CONTAINS              | NUM       | LN 11627-7<br>Amniotic Fluid<br>Index     | UCUM cm<br>centimeter             | 300        |            | AFI                                  |
|        |    |                       |           | LN 11624-4 First<br>Quadrant<br>Diameter  |                                   |            | 12008      | Q1                                   |
|        |    |                       |           | LN 11626-9<br>Second Quadrant<br>Diameter |                                   |            |            | Q2                                   |
|        |    |                       |           | LN 11625-1 Third<br>Quadrant<br>Diameter  |                                   |            |            | Q3                                   |

|  |  |  |  |                 |  |  |  |    |
|--|--|--|--|-----------------|--|--|--|----|
|  |  |  |  | LN 11623-6      |  |  |  |    |
|  |  |  |  | Fourth Quadrant |  |  |  | Q4 |
|  |  |  |  | Diameter        |  |  |  |    |

**9.1.1.12 Pelvis and Uterus Section (TID 5015)**




**Figure 9.1-11**  
**TEMPLATE HIERARCHY OF PELVIS AND UTERUS SECTION IN OB-GYN SR**

**Table 9.1-12**  
**PELVIS AND UTERUS SECTION IN OB-GYN SR**

|          | NL  | REL      | VT        | Concept Name                    | Unit / CODE Value     | Ref TID | Ref CID | Comments  |
|----------|-----|----------|-----------|---------------------------------|-----------------------|---------|---------|---|
| A-18     | >   | CONTAINS | CONTAINER | DCM 125011<br>Pelvis and Uterus |                       | 5015    |         |   |
| A-18-1   | >>  | CONTAINS | CONTAINER | SRT T-83000<br>Uterus           |                       | 5016    |         | TID 5016(LWH Volume Group) is included. Group Name is "Uterus". |
| A-18-1-1 | >>> | CONTAINS | NUM       | LN 11865-3<br>Uterus Width      | UCUM cm<br>centimeter | 300     |         | This row is inserted as part of TID 300 invocation.             |
|          |     |          |           | LN 11842-2<br>Uterus Length     |                       |         |         | This row is inserted as part of TID 300 invocation.             |

|            |      |                       |      |   |  |       |  |
|------------|------|-----------------------|------|---|--|-------|--|
|            |      |                       |      | LN 11859-6<br>Uterus Height   |  |       | This row is inserted as part of TID 300 invocation.  |
| A-18-1-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation  | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value"  |
| A-18-1-2   | >>>  | CONTAINS              | NUM  | LN 33192-6<br>Uterus Volume   | UCUM cm3<br>Cubic<br>centimeter                  |       | This row is inserted as part of TID 300 invocation. X8 automatically calculates the volume based on Length, Width and Height measurements.                           |
| A-18-2     | >>   | CONTAINS              | NUM  | LN 11961-0<br>Cervix Length<br><br>LN 12145-9<br>Endometrium<br>Thickness | UCUM cm<br>centimeter                            | 12011 | This measurement is from CID 12011. This row is inserted as part of TID 300 invocation. Only Cervix Length and Endometrium Thickness from CID 12011 will be present. |
| A-18-2-1   | >>>  | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation  | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value"  |

9.1.1.13 Ovaries Section (TID 5012)



**Figure 9.1-12**  
**TEMPLATE HIERARCHY OF OVARIES SECTION IN OB-GYN SR**


**Table 9.1-13**  
**OVARIES SECTION IN OB-GYN SR**

|          | NL  | REL             | VT        | Concept Name                | Unit / CODE Value   | Ref TID | Ref CID | Comments  |
|----------|-----|-----------------|-----------|-----------------------------|---------------------|---------|---------|---|
| A-19     | >   | CONTAINS        | CONTAINER | DCM 121070 Findings         |                     | 5012    |         |   |
| A-19-1   | >>  | HAS CONCEPT MOD | CODE      | SRT G-C0E3 Finding Site     | SRT T-87000 "Ovary" |         |         | DT (T-87000, SRT, "Ovary")  |
| A-19-2   | >>  | CONTAINS        | CONTAINER | SRT T-87000 Ovary           |                     | 5016    |         | TID 5016(LWH Volume Group) is included. Left ovary volume, length and width measurements are inserted. Group Name is "Ovary". |
| A-19-2-1 | >>> | CONTAINS        | NUM       | LN 11829-9 Left Ovary Width | UCUM cm centimeter  | 300     |         | This row is inserted as part of TID 300 invocation.   |

|            |      |                       |           |                                  |  |       |  |
|------------|------|-----------------------|-----------|----------------------------------|--|-------|--|
|            |      |                       |           | LN 11840-6 Left<br>Ovary Length  |  |       | This row is inserted as part of TID 300 invocation.  |
|            |      |                       |           | LN 11857-0 Left<br>Ovary Height  |  |       | This row is inserted as part of TID 300 invocation.  |
| A-19-2-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE      | DCM 121401<br>Derivation         | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 | If user selects Avg. from the Report, this value will be a "Mean". Else this value will be a "Best value"                                  |
| A-19-2-2   | >>>  | CONTAINS              | NUM       | LN 12164-0 Left<br>Ovary Volume  | UCUM cm3<br>Cubic<br>centimeter                  |       | This row is inserted as part of TID 300 invocation. X8 automatically calculated the volume based on Length, Width and Height measurements. |
| A-19-3     | >>   | CONTAINS              | CONTAINER | SRT T-87000<br>Ovary             |  | 5016  | Similarly TID 5016(LWH Volume Group) is included for Right Ovary Volume, length and width measurements.                                    |
| A-19-3-1   | >>>  | CONTAINS              |           | LN 11830-7 Right<br>Ovary Width  | UCUM cm<br>centimeter                            |       |  |
|            |      |                       |           | LN 11841-4 Right<br>Ovary Length |  |       |  |
|            |      |                       | NUM       | LN 11858-8 Right<br>Ovary Height |  |       |  |
| A-19-3-1-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE      | DCM 121401<br>Derivation         | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 |  |
| A-19-3-2   | >>>  | CONTAINS              | NUM       | LN 12165-7 Right<br>Ovary Volume | UCUM cm3<br>Cubic                                |       |  |

|  |  |  |  |  |            |  |  |  |
|--|--|--|--|--|------------|--|--|--|
|  |  |  |  |  | centimeter |  |  |  |
|--|--|--|--|--|------------|--|--|--|

**9.1.1.14 Follicles Section - Left (TID 5013)**




**Figure 9.1-13**  
**TEMPLATE HIERARCHY OF FOLLICLES SECTION (LEFT) IN OB-GYN SR**

**Table 9.1-14**  
**FOLLICLES SECTION (LEFT) IN OB-GYN SR**

|        | NL | REL                   | VT        | Concept Name  | Unit / CODE<br>Value                 | Ref<br>TID | Ref<br>CID | Comments                                 |
|--------|----|-----------------------|-----------|---|--------------------------------------|------------|------------|--|
| A-20   | >  | CONTAINS              | CONTAINER | DCM 121070<br>Findings                                |                                      | 5013       |            |  |
| A-20-1 | >> | HAS<br>CONCEPT<br>MOD | CODE      | SRT G-C0E3<br>Finding Site                            | SRT T-87600<br>"Ovarian<br>Follicle" |            |            | DT (T-87600, SRT,<br>"Ovarian Follicle") |
| A-20-2 | >> | HAS<br>CONCEPT<br>MOD | CODE      | SRT G-C171<br>Laterality                              | SRT G-A101<br>Left                   |            |            | EV (G-A101, SRT,<br>"Left")              |
| A-20-3 | >> | CONTAINS              | NUM       | LN 11879-4<br>Number of<br>follicles in left<br>ovary | UCUM 1 no<br>units                   |            |            |  |
| A-20-4 | >> | CONTAINS              | CONTAINER | DCM 125007<br>Measurement                             |                                      | 5014       |            | TID 5014 (Follicle<br>Measurement)       |

|                |      |                       |      |                                 |  |       |   |
|----------------|------|-----------------------|------|---------------------------------|--|-------|---|
|                |      |                       |      | Group                           |  |       | Group) is included.   |
| A-20-4-1       | >>>  | HAS OBS<br>CONTEXT    | TEXT | DCM 12510<br>Identifier         |  |       | Uses numbers "1,<br>"2, "3 . up to "12 to<br>identify the follicle.   |
| A-20-4-2       | >>>  | CONTAINS              | NUM  | SRT GD705<br>Volume             | UCUM cm3<br>Cubic<br>centimeter                  |       | This is inserted as<br>part of TID 300<br>invocation. X8<br>automatically<br>calculates the<br>volume based on<br>the follicle diameter |
| A-20-4-3       | >>>  | CONTAINS              | NUM  | LN 11793-7<br>Follicle Diameter | UCUM cm<br>centimeter                            | 300   | This is inserted as<br>part of TID 300<br>invocation. This<br>name for X8 is<br>[1],[2],[3],...[12].                                    |
| A-20-4-3-<br>1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation        | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 | If user selects Avg.<br>from the Report,<br>this value will be a<br>"Mean". Else this<br>value will be a<br>"Best value"                |

### 9.1.1.15Follicles Section – Right (TID 5013)



**Figure 9.1-14**  
**TEMPLATE HIERARCHY OF FOLLICLES SECTION (RIGHT) IN OB-GYN SR**

**Table 9.1-15**  
**FOLLICLES SECTION (RIGHT) IN OB-GYN SR**

|          | NL  | REL                   | VT        | Concept Name   | Unit / CODE<br>Value                 | Ref<br>TID | Ref<br>CID | Comments   |
|----------|-----|-----------------------|-----------|--|--------------------------------------|------------|------------|--|
| A-21     | >   | CONTAINS              | CONTAINER | DCM 121070<br>Findings                                 |                                      | 5013       |            |  |
| A-21-1   | >>  | HAS<br>CONCEPT<br>MOD | CODE      | SRT G-C0E3<br>Finding Site                             | SRT T-87600<br>"Ovarian<br>Follicle" |            |            | DT (T-87600, SRT,<br>"Ovarian Follicle")   |
| A-21-2   | >>  | HAS<br>CONCEPT<br>MOD | CODE      | SRT G-C171<br>Laterality                               | SRT G-A100<br>"Right"                |            |            | EV (G-A100, SRT,<br>"Right")   |
| A-21-3   | >>  | CONTAINS              | NUM       | LN 11880-2<br>Number of<br>follicles in right<br>ovary | UCUM 1 no<br>units                   |            |            |  |
| A-21-4   | >>  | CONTAINS              | CONTAINER | DCM 125007<br>Measurement<br>Group                     |                                      | 5014       |            | Similar TID<br>5014(Follicle<br>Measurement<br>Group) is included<br>for follicles in right<br>ovary diameter and<br>volume. |
| A-21-4-1 | >>> | HAS OBS<br>CONTEXT    | TEXT      | DCM 12510<br>Identifier                                |                                      |            |            |  |
| A-21-4-2 | >>> | CONTAINS              | NUM       | SRT G-D705<br>Volume                                   | UCUM cm3<br>Cubic<br>centimeter      | 300        |            |  |

|            |      |                       |      |                                 |  |       |  |
|------------|------|-----------------------|------|---------------------------------|--|-------|--|
| A-21-4-3   | >>>  | CONTAINS              | NUM  | LN 11793-7<br>Follicle Diameter | UCUM cm<br>centimeter                            |       |  |
| A-21-4-3-1 | >>>> | HAS<br>CONCEPT<br>MOD | CODE | DCM 121401<br>Derivation        | SRT R-002E1<br>Best value<br>SRT R-00317<br>Mean | /3627 |  |

## 9.2 DCMR Context Groups used in SONOACE X8

### 9.2.1 Gestational Age Equations and Tables (Context Group 12013)

**Table 9.2-1  
GESTATIONAL AGE EQUATIONS AND TABLES**

| <b>Coding Scheme Designator (0008,0102)</b> | <b>Code Value (0008,0100)</b> | <b>Code Meaning (0008,0104)</b> |
|---|-------------------------------|---------------------------------|
| LN  | 11889-3                       | AC, Campbell 1975               |
| LN  | 11892-7                       | AC, Hadlock 1984                |
| LN  | 33076-1                       | AC, Shinozuka 1996              |
| LN  | 11902-4                       | BPD, Hadlock 1984               |
| LN  | 33538-0                       | BPD, Hansmann 1986              |
| LN  | 11905-7                       | BPD, Jeanty 1984                |
| LN  | 11906-5                       | BPD, Kurtz 1980                 |
| LN  | 33082-9                       | BPD, Osaka 1989                 |
| LN  | 11907-3                       | BPD, Sabbagha 1978              |
| LN  | 33084-5                       | BPD, Shinozuka 1996             |
| LN  | 33086-0                       | BPD-oi, Chitty 1997             |
| LN  | 33087-8                       | BPD-oo, Chitty 1997             |
| LN  | 33088-6                       | Clavical length, Yarkoni 1985   |
| LN  | 11910-7                       | CRL, Hadlock 1992               |
| LN  | 33540-6                       | CRL, Hansmann 1986              |
| LN  | 11913-1                       | CRL, Nelson 1981                |
| LN  | 33093-6                       | CRL, Osaka 1989                 |
| LN  | 33094-4                       | CRL, Rempen 1991                |

|    |         |                            |
|----|---------|----------------------------|
| LN | 11914-9 | CRL, Robinson 1975         |
| LN | 33095-1 | CRL, Shinozuka 1996        |
| LN | 33098-5 | FL, Chitty 1997            |
| LN | 11920-6 | FL, Hadlock 1984           |
| LN | 33541-4 | FL, Hansmann 1986          |
| LN | 11922-2 | FL, Hohler 1982            |
| LN | 11923-0 | FL, Jeanty 1984            |
| LN | 33101-7 | FL, Osaka 1989             |
| LN | 33102-5 | FL, Shinozuka 1996         |
| LN | 11928-9 | GS, Hellman 1969           |
| LN | 33107-4 | GS, Nyberg 1992            |
| LN | 33108-2 | GS, Tokyo 1986             |
| LN | 33110-8 | HC measured, Chitty 1997   |
| LN | 33111-6 | HC derived, Chitty 1997    |
| LN | 11932-1 | HC, Hadlock 1984           |
| LN | 33543-0 | HC, Hansmann 1986          |
| LN | 11936-2 | Humerus, Jeanty 1984       |
| LN | 33117-3 | Humerus Length, Osaka 1989 |
| LN | 33120-7 | OFD, Hansmann 1986         |
| LN | 11941-2 | Tibia, Jeanty 1984         |
| LN | 11944-6 | Ulna, Jeanty 1984          |

## 9.2.2 OB Fetal Body Weight Equations and Tables (Context ID 12014)

**Table 9.2-2**

### **OB FETAL BODY WEIGHT EQUATIONS AND TABLES**

| <b>Coding Scheme Designator (0008,0102)</b> | <b>Code Value (0008,0100)</b> | <b>Code Meaning (0008,0104)</b>      |
|---|-------------------------------|--------------------------------------|
| LN  | 11756-4                       | EFW by AC, Campbell 1975             |
| LN  | 11738-2                       | EFW by AC, BPD, Hadlock 1984         |
| LN  | 11735-8                       | EFW by AC, BPD, FL, Hadlock 1985     |
| LN  | 11732-5                       | EFW by AC, BPD, FL, HC, Hadlock 1985 |

|    |         |                                 |
|----|---------|---------------------------------|
| LN | 11751-5 | EFW by AC, FL, Hadlock 1985     |
| LN | 11746-5 | EFW by AC, FL, HC, Hadlock 1985 |
| LN | 33139-7 | EFW by BPD, TTD, Hansmann 1986  |
| LN | 11739-0 | EFW by AC and BPD, Shepard 1982 |
| LN | 33140-5 | EFW by BPD, FTA, FL, Osaka 1990 |

### 9.2.3 Fetal Growth Equations and Tables (Context ID 12015)

**Table 9.2-3  
FETAL GROWTH EQUATIONS AND TABLES**

| <b>Coding Scheme Designator (0008,0102)</b> | <b>Code Value (0008,0100)</b> | <b>Code Meaning (0008,0104)</b>    |
|---|-------------------------------|------------------------------------|
| LN  | 33145-4                       | AC by GA, ASUM 2000                |
| LN  | 33146-2                       | AC by GA, Hadlock 1984             |
| LN  | 33147-0                       | AC (measured) by GA, Chitty 1994   |
| LN  | 33546-3                       | AC (derived) by GA, Chitty 1994    |
| LN  | 33149-6                       | AC by GA, Shinozuka 1996           |
| LN  | 33151-2                       | BPD by GA, ASUM 2000               |
| LN  | 33198-3                       | BPD by GA, Hadlock 1984            |
| LN  | 33556-2                       | BPD outer-inner by GA, Chitty 1994 |
| LN  | 33152-0                       | BPD outer-outer by GA, Chitty 1994 |
| LN  | 33156-1                       | BPD by GA, Shinozuka 1996          |
| LN  | 33161-1                       | CRL by GA, Shinozuka 1996          |
| LN  | 33164-5                       | Fibula by GA, Jeanty 1983          |
| LN  | 33165-2                       | FL by GA, ASUM 2000                |
| LN  | 33166-0                       | FL by GA, Hadlock 1984             |
| LN  | 33167-8                       | FL by GA, Chitty 1994              |
| LN  | 33170-2                       | FL by GA, Shinozuka 1996           |
| LN  | 33172-8                       | HC by GA, ASUM 2000                |
| LN  | 33173-6                       | HC by GA, Hadlock 1984             |
| LN  | 33174-4                       | HC derived by GA, Chitty 1994      |
| LN  | 33177-7                       | Humerus Length by GA, ASUM 2000    |
| LN  | 33178-5                       | OFD by GA, ASUM 2000               |

|    |         |                           |
|----|---------|---------------------------|
| LN | 33180-1 | Radius by GA, Jeanty 1983 |
| LN | 33181-9 | TCD by GA Goldstein 1987  |

#### 9.2.4 Estimated Fetal Weight Percentile Equations and Tables (Context ID 12016)

**Table 9.2-4**

**ESTIMATED FETAL WEIGHT PERCENTILE EQUATIONS AND TABLES**

| <b>Coding Scheme<br/>Designator<br/>(0008,0102)</b> | <b>Code Value<br/>(0008,0100)</b> | <b>Code Meaning<br/>(0008,0104)</b> |
|---|-----------------------------------|-------------------------------------|
| LN  | 33183-5                           | FWP by GA, Hadlock 1991             |
| LN  | 33184-3                           | FWP by GA, Williams, 1982           |
| LN  | 33189-2                           | FWP by GA, Brenner 1976             |